Technical Data Sheet: Raji-Luc2

ATCC® Number	CCL-86-LUC2™
Organism	Homo sapiens
Tissue/Disease Source	Burkitt's lymphoma
Product Description	This luciferase expressing cell line was derived from Raji cell line by transduction with lentiviral vector encoding firefly luciferase gene (luc2) and subsequently through single cell cloning. • Signal noise ratio: ≥ 1,000 • Naturally expresses high levels of CD19 (verified at ATCC) • Raji-Luc2 has been used as a target cancer cell for <i>in vitro</i> killing assay by CD19 CAR-T cells (tested at ATCC) and is expected to also work for CD20 CAR-T cells.
Application	Excellent signal/background ratio and stable luciferase expression make this cell line ideal for in vitro study of CD19 specific CAR-T cells. It also can be used in cell-based assays for cancer research.

In vitro CAR-T killing cancer bioluminescence assay




Figure 1: CD19 CAR-T *in vitro* **killing assay of Raji-Luc2 measured using luminescence.** Raji-Luc2 cells (5 x 10³) were seeded into a 96-well plate and were used as target cells for either CD19 CAR-T or Mock CAR-T (control) from the same donor which were seeded at various ratios of CAR-T cells to target Raji-Luc2 cells (1:1, 2:1, 5:1, and 10:1). After 24 hours of co-culture, Bright-Glo (Promega) was added to the indicated wells. The luminescence of the plate was read within 10 minutes using a luminescence plate reader and determined to have a dose-dependent specific killing with CD19 CAR-T cells which was greater than the non-specific killing observed with mock CAR-T cells. (* = significant difference, ns = not significant using unpaired t test, with a single pooled variance).

Cell Morphology




Figure 2: Cell morphology of Raji parental and Raji-Luc2. Cells were maintained in ATCC recommended culture conditions. Cell morphology was observed under microscopy and images were captured by digital camera.

Luciferase Expression




Figure 3: Linearity of luminescence of Raji-Luc2. Cells were seeded in a 96-well plate at indicated cell numbers per well, and Bright-Glo (Promega) was added to the indicated wells. The luminescence of the plate was read within 10 minutes using a luminescence plate reader and was determined to have a linear correlation of bioluminescence intensity with cell numbers.

© 2022 American Type Culture Collection. The ATCC trademark and trade name, and any other trademarks listed in this publication are owned by the American Type Culture Collection unless indicated otherwise.

These products are for laboratory use only. Not for human or diagnostic use. ATCC products may not be resold, modified for resale, used to provide commercial services, or to manufacture commercial products without prior ATCC written approval.