

## **ATCC Medium: 37 Chromatium Medium**

### **Solution 1**

CaCl <sub>2</sub> .....	0.8 g
DI Water.....	1000 ml

Adjust the pH to 6.6-6.8. The medium should be dispensed into sterile 16x125 screw cap tubes. Fill the tubes to about 1/8 of an inch from the top.

Autoclave at 121°C

**Note:** Media may be dispensed aseptically to avoid boil over. If not, I would recommend using rubber screw caps.

### **Solution 2**

Heavy Metal Solution (see below).....	50.0 ml
Vitamin B <sub>12</sub> Solution (see below).....	3.0 ml
Vitamin Solution (see below).....	15.0 ml
KH <sub>2</sub> PO <sub>4</sub> .....	1.0 g
KCl.....	1.0 g
NH <sub>4</sub> Cl.....	0.8 g
MgCl <sub>2</sub> x 6H <sub>2</sub> O.....	0.8 g
Sodium Ascorbate.....	2.4 g
DI Water.....	32.0 ml

### **Solution 3**

NaHCO <sub>3</sub> .....	4.5 g
DI Water.....	900.0 ml

Combine Solutions 2 and 3. Adjust the pH to 6.6-6.8. The medium should be dispensed into sterile 16x125 screw cap tubes. Fill the tubes to about 1/8 of an inch from the top.

Autoclave @ 121C.

**Note:** Media may be dispensed aseptically to avoid boil over. If not, I would recommend using rubber screw caps.

### **Solution 4**

**This may be added by the bacteriologist.**

Na <sub>2</sub> S x 9H <sub>2</sub> O.....	3.0 g
DI Water.....	200.0 ml

Gas under N<sub>2</sub>. Dispense into Hungate tubes or bottles. Generally 4-5 test tubes (8-10ml each is sufficient). Autoclave @ 121C

### **Heavy Metal Solution**

Ethylene Diamine Tetraacetate.....	1.5 g
------------------------------------	-------

Modified Hogland trace elements solution (see below).....	6.0 ml
FeSO <sub>4</sub> x 7H <sub>2</sub> O.....	0.2 g
ZnSO <sub>4</sub> x 7H <sub>2</sub> O.....	0.1 g
MnCl <sub>2</sub> x 4H <sub>2</sub> O.....	0.02 g
DI Water.....	1000 ml

\*Dissolve EDTA first

**Vitamin B<sub>12</sub> Solution**

Vitamin B <sub>12</sub> .....	2.0 mg
DI Water.....	100.0 ml

**Vitamin Solution**

**(May substitute ATCC Vitamin Supplement)**

Biotin.....	0.2 mg
Nicotinic Acid.....	2.0mg
Thiamin.....	1.0 mg
p-Aminobenzoic acid.....	0.1mg
Pantothenic acid.....	0.5 mg
Pyrodoxamine HCl.....	5.0 mg
DI Water.....	100.0 ml

**Modified Hogland Trace Element Solution**

**(May substitute ATCC Trace Minerals)**

AlCl <sub>3</sub> .....	1.0 g
KI.....	1.0 g
KBr.....	0.5 g
LiCl.....	0.5 g
MnCl <sub>2</sub> x 4H <sub>2</sub> O.....	7.0 g
H <sub>3</sub> BO <sub>3</sub> .....	11.0 g
ZnCl <sub>2</sub> .....	1.0 g
CuCl <sub>2</sub> .....	1.0 g
NiCl <sub>2</sub> .....	1.0 g
CoCl <sub>2</sub> .....	1.0 g
SnCl <sub>2</sub> x 2H <sub>2</sub> O.....	0.5 g
BaCl <sub>2</sub> .....	0.5 g
Na <sub>2</sub> MoO <sub>4</sub> .....	0.5 g
NaVO <sub>3</sub> x H <sub>2</sub> O.....	0.1 g
Se-Salt (Na selenite).....	0.5 g
DI Water.....	3600 ml

Each salt is dissolved in distilled water before mixing; pH of each solution is adjusted to below 7.0. The final pH is adjusted to 3-4. Before use, the solution is thoroughly mixed. The total final volume is 3.6 liters. The pH of the final solution is adjusted to pH 3-4. The flaky yellow precipitate which is formed after mixing transforms after a few days into a very fine white precipitate. Before use the

solution is mixed thoroughly.