

## Protocol for OAT1 HEK 293T/17 Uptake of 5-Carboxyfluorescein (5-CF)

There is an abundance of OAT1 transporter in kidney tissue. However primary cells lose OAT1 expression after just a few days in culture. To facilitate hOAT1 *in-vitro* toxicity studies we have generated OAT1 HEK 293T/17 (ATCC® CRL-11268G-1™) cell line by stably expressing hOAT1 into the HEK 293T/17 (ATCC® CRL-11268™) cell line.

Functionality of the OAT1 HEK 293T/17 was demonstrated by uptake of 5-Carboxyfluorescein (5-CF). 5-CF is a fluorescent in the FITC channel and is readably taken up by the OAT1 transporter. After uptake, the cells can be visualized under a fluorescent microscope or lysed and read on a fluorescent plate reader.

### Materials:

Material	Company	Catalogue Number
OAT1 HEK 293T/17	ATCC	CRL-11268G-1
HEK 293T/17	ATCC	CRL-11268
HBSS	Corning	21-023-CV
5-CF	Sigma	86826
Black 96-well plates	Corning	354640
M-Per lysis buffer	Thermo Fisher	78501
DMEM	ATCC	30-2002
FBS	ATCC	30-2020


### Protocol:

#### A. Seeding cells

1. Thaw or split both HEK 293T/17 and OAT1 HEK 293T/17.
2. Resuspend the each cell line at a concentration of  $5 \times 10^5$  cells/mL in DMEM with 10% FBS.
3. Plate 200 $\mu$ L of cells into each well of a 96-well poly-D-lysine coated plate ensuring that both cells have a minimum of three replicates each.
4. Incubate the cells for 18 – 28 hours at 37°C and 5% CO<sub>2</sub> incubator.

#### B. 5-CF assay

1. Make sure the cells are more than 90% confluent before you start the assay.
2. Wash the cells three times with 200  $\mu$ L warm (37°C) HBSS<sup>1</sup>.
3. Incubate the third wash at 37°C and 5% CO<sub>2</sub> for 10 minutes.
4. Add 200  $\mu$ L of 150  $\mu$ M of 5-CF in warm HBSS and incubate at 37°C and 5% CO<sub>2</sub> for 20 minutes<sup>2</sup>.
5. Wash three times with cold (4°C) HBSS.
6. (Optional) The last wash can remain on the cells for up to 10 minutes to visualize the cells in the FITC channel.
7. Remove the wash and add 100  $\mu$ L of M-Per lysis buffer.
8. Incubate for 5 minutes at room temperature protected from light.
9. Read on a fluorescent plate reader at 490<sub>ex</sub> and 510-580<sub>em</sub>.
10. Calculate the signal of OAT over parental to determine the relative ratio. Or the parental can be used to subtract background.




**Figure1. Uptake of 5-CF by OAT1 HEK 293T/17 (ATCC<sup>®</sup> CRL-11268G-1<sup>™</sup>)**

OAT1 HEK 293T/17 cells after taking up 5-CF visualized using a Nikon fluorescent microscope.

<sup>1</sup> HEK 293T/17 cells are loosely adherent. When adding or removing liquid try not to disturb the cells.

<sup>2</sup> 5-CF is light sensitive and care should be taken with the compound and once the cells take up the compound to protect from light.



**Figure1. Graph of OAT1 HEK 293T/17 absolute 5-CF signal over three passages.**

OAT1 HEK 293T/17 and HEK 293T/17 absolute signal of 5-CF after lysis and monitored on a Promega GlowMax over three different passages.