





Product Sheet

# Chondrocyte Differentiation Tool (ATCC® PCS-500-051™)

Please read this FIRST

 Storage Temp.  
**-20°C (or -70°C for long-term storage)**

---

 Biosafety Level  
**1**

## Description

### Product Description:

The Chondrocyte Differentiation Tool (ATCC® PCS-500-051) is a complete differentiation medium designed to induce chondrogenesis in actively proliferating Adipose-Derived Mesenchymal Stem Cells (ATCC® PCS-500-011) with high efficiency.

**Formulation:** Proprietary

**Volume:** 100 mL

## Components

**Formulation:** Proprietary

**Volume:** 100 mL\*

\*The Osteocyte Differentiation Tool provides enough medium for differentiation of ~ 1 million cells when plated at a recommended density of 18,000 viable cells/cm<sup>2</sup> in a 6-well tissue culture format.

## Directions for Use

Antimicrobials and phenol red are not required but may be added to the Chondrocyte Differentiation Tool if desired prior to use. The recommended volume of each **optional** component to be added to Chondrocyte Differentiation Tool is summarized in Table 1.

**Table 1.** Optional Addition of Antimicrobials/Antimycotics and Phenol Red per 100 mL of Medium

Component	Volume	Final Concentration
Gentamicin-Amphotericin B Solution	0.1 mL	Gentamicin: 10 µg/mL Amphotericin B: 0.25 µg/mL
Penicillin-Streptomycin-Amphotericin B Solution	0.1 mL	Penicillin: 10 Units/mL Streptomycin: 10 µg/mL Amphotericin B: 25 ng/mL
Phenol Red	0.1 mL	33 µM

### Additional Materials Needed for Alginate Encapsulation (not provided):

1. Sodium alginate (e.g., Sigma-Aldrich #71238)
2. 3 mL syringe
3. 27-gauge needle
4. Steriflip® Filter Unit, Millipore Corporation (Vacuum-driven 50 mL filtration system)
5. Wide bore pipette tip
6. 150 mM NaCl solution
7. 100 mM CaCl<sub>2</sub> solution (sterile)
8. Small magnetic stir bar (sterile)
9. 250 mL Beaker (sterile)
10. Forceps or stir bar extractor (sterile)
11. Magnetic stir plate

### Chondrocyte Differentiation

Chondrocyte differentiation requires that the cells be grown in a three-dimensional aggregate cell culture. Micromass culture can be used; but, for the best results, ATCC recommends the use of a matrix, such as alginate, to provide a scaffold for the deposition of proteoglycans. The following procedure demonstrates differentiation of ~2.7 x 10<sup>7</sup> cells seeded to four wells of a 48-well tissue culture plate, using 25 mL of the Chondrocyte Differentiation Tool and alginate encapsulation.

### Preparation of a 1.5% (w/v) Alginate Solution

1. Add 0.15 g of alginate to 10 mL of 150 mM NaCl while stirring rapidly or vortexing to minimize clumping.
2. Agitate the solution on a rocker for at least 2 hr and up to 16 hr at room temperature to completely solubilize the alginate.
3. Sterilize the solution using a 0.22 µm filter. Store at 4°C to 8°C for up to week.

### Preparation of Differentiation Medium

1. Pre-warm the Chondrocyte Differentiation Tool to 37°C in a water bath.

American Type Culture Collection  
PO Box 1549  
Manassas, VA 20108 USA  
[www.atcc.org](http://www.atcc.org)

800.638.6597 or 703.365.2700  
Fax: 703.365.2750  
Email: [Tech@atcc.org](mailto:Tech@atcc.org)


Or contact your local distributor



Product Sheet


## Chondrocyte Differentiation Tool (ATCC® PCS-500-051™)

Please read this FIRST



Storage Temp.  
**-20°C (or -70°C for long-term storage)**

---



Biosafety Level  
**1**

2. Once thawed, store the remaining Chondrocyte Differentiation Tool in the dark at 2°C-8°C for later use. When stored under these conditions, the differentiation media is stable for up to three weeks.

**Note:** For procedures using less than the full 100 mL volume of the Chondrocyte Differentiation Tool, the medium can be dispensed into appropriate aliquots (e.g., for use with this protocol, we recommend dispensing into 25 mL aliquots). Unused portions can be refrozen once without loss of efficiency; however, multiple freeze/thaw cycles are NOT recommended.

### Preparing Cells for Chondrocyte Differentiation

1. Follow the instructions for the growth of Adipose-Derived Mesenchymal Stem Cells (ATCC® PCS-500-011) using Mesenchymal Stem Cell Basal Medium (ATCC® PCS-500-030) supplemented with Mesenchymal Stem Cell Growth Kit – Low Serum (ATCC® PCS-500-040) components.
2. Expand the cells as needed for experimental design, but do not passage more than four (4) times prior to initiating chondrocyte differentiation.

**Note:** It will take approximately 10 to 15 T-75 flasks at 70%-80% confluence to obtain  $2.7 \times 10^7$  cells, which is enough to seed four wells of a 48-well tissue culture plate using the alginate encapsulation method described below.

1. Cells should be collected and counted when the culture is 70%-80% confluent and actively proliferating.
2. After centrifugation at 150 x g for 3-5 minutes, discard the supernatant and resuspend the cell pellet ( $2.5 \times 10^7$  cells) in 800  $\mu$ L of the 1.5% (w/v) Alginate Solution. The final volume should be ~1.0 mL.

**Note:** The Alginate Solution must not be diluted lower than 1.2% (w/v) by the addition of the cells. If you have fewer cells, adjust the volume of the prepared 1.5% Alginate Solution used in order to maintain this ratio. Likewise, if you have more cells, scale proportionally.

1. Gently mix the cell-alginate suspension by pipetting up and down; taking care not to introduce bubbles into the solution.
2. Proceed to Step 1 (below) for alginate encapsulation of cells.

### Alginate Encapsulation of Cells: Formation of Chondrogenic Microbeads

1. Transfer 75 mL of a 100 mM sterile CaCl<sub>2</sub> solution to a sterile 250 mL beaker containing a sterile stir bar.
2. Create a gentle funnel in the CaCl<sub>2</sub> solution at room temperature on a stir plate.
3. Transfer the alginate-cell suspension with a 3 mL syringe fitted with a 27 gauge needle.
4. Rapidly dispense the alginate-cell suspension into the CaCl<sub>2</sub> solution to form chondrogenic microbeads.
5. Allow the chondrogenic microbeads to stir for 10 minutes to solidify (cure) the alginate.
6. Remove the beaker from the stir plate and allow the chondrogenic microbeads to settle on the bottom.
7. Transfer the chondrogenic microbead solution into a 50 mL conical tube and attach to the vacuum-driven Steriflip® Filter Unit. Immediately break the vacuum as soon as the liquid is removed to prevent damage to the beads.
8. Resuspend the chondrogenic microbeads in 2 mL of pre-warmed Chondrocyte Differentiation Tool.
9. Aseptically transfer enough chondrogenic microbeads to cover the bottom surfaces of four wells of a 48-well plate (about 0.5 mL per well). This set up will yield  $\sim 6.75 \times 10^6$  encapsulated cells/well.
10. Allow the chondrogenic microbeads to settle to the bottom of the wells. To remove residual CaCl<sub>2</sub>, wash the chondrogenic microbeads by replacing the medium in each well twice with 0.5 mL Chondrocyte Differentiation Tool. Add 0.5 mL Chondrocyte Differentiation Tool to each well after the last wash.
11. Incubate the cells at 37°C with 5% CO<sub>2</sub> for 2-3 days before renewing the medium.
12. When ready to renew the medium, retrieve the Chondrocyte Differentiation Tool from storage and transfer the required volume to a sterile tube. (For 4 wells in a 48-well plate, this volume would be 2 mL).
13. Warm the aliquot of Chondrocyte Differentiation Tool to 37°C in a water bath.
14. Carefully remove the spent medium, taking great care not to disturb or aspirate the chondrogenic microbeads.
15. Add 0.5 mL of fresh, pre-warmed Chondrocyte Differentiation Tool to each well.
16. Incubate the cells at 37°C with 5% CO<sub>2</sub> for 2-3 days before renewing the medium.
17. Repeat steps 12 through 16 every 2-3 days until the cells have been exposed to the Chondrocyte Differentiation Tool for a total of 21 days.
18. Cells can be used at any phase of chondrocyte differentiation as predicated upon experimental design. To confirm calcium accumulation, cells can be fixed and stained with Alcian Blue (not provided).
- 19.

American Type Culture Collection  
PO Box 1549  
Manassas, VA 20108 USA  
[www.atcc.org](http://www.atcc.org)

800.638.6597 or 703.365.2700  
Fax: 703.365.2750  
Email: [Tech@atcc.org](mailto:Tech@atcc.org)

Or contact your local distributor

### ATCC Warranty


The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the



Product Sheet


## Chondrocyte Differentiation Tool (ATCC® PCS-500-051™)

Please read this FIRST



Storage Temp.  
**-20°C (or -70°C for  
long-term  
storage)**

---



Biosafety Level  
**1**

product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

---

### Disclaimers

---

This product is intended for laboratory research purposes only. It is not intended for use in humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at [www.atcc.org](http://www.atcc.org)

Additional information on this culture is available on the ATCC web site at [www.atcc.org](http://www.atcc.org).

© ATCC 2016. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [08/08]

American Type Culture Collection  
PO Box 1549  
Manassas, VA 20108 USA  
[www.atcc.org](http://www.atcc.org)

800.638.6597 or 703.365.2700  
Fax: 703.365.2750  
Email: [Tech@atcc.org](mailto:Tech@atcc.org)

Or contact your local distributor