



The Battle Against Antimicrobial Resistance: ATCC® is your Resource for AMR strains from Drug Discovery to Clinical testing

Victoria Knight-Connoni, PhD
Head of Content Development, ATCC

Joe Petrone, PhD
Senior Bioinformatician, ATCC



About Us



ATCC is a global leader in providing authenticated, high-quality biological resources and standards for industry, academia, and government.

Founded in 1925, ATCC is a private, nonprofit, global biological resource center and standards organization that provides scientists with the biomaterials and resources they need to conduct critical life science research.



Outline

- The global burden of antimicrobial resistance
- The rise and spread of resistance
- Challenges in antibacterial drug discovery
- ATCC solutions
 - Over 200 multidrug-resistance clinical isolates with source metadata and AST data
 - Complete *de novo* hybrid genome assemblies with annotated antibiotic resistance genes



The Global Burden

- **Global Impact:** At least 1.14 million deaths were directly caused by AMR in 2021, with nearly 5 million associated deaths globally (*Lancet 2024*).
- **U.S. Burden:** Over 2.8 million AMR infections and 48,000 deaths annually (*CDC 2022*).
- **Economic Cost:** Treating the six most common AMR infections in U.S. healthcare costs over \$4.6 billion per year (*Nelson et al 2025 Clin Infect Dis 72: S17-S26*).

Factors Increasing Susceptibility to Infection and AMR




Age

- Very young
- Older adults (>65 years)



Illness

- Underlying chronic health conditions
- Compromised immune system
- Recent Infection



Nutrition

- Malnutrition
- Obesity




Immuno-suppressive drugs or treatments

- Chemotherapy, radiation, and immunosuppressants
- Long-term antibiotic use

The Rise of Resistance

Significant events in antibiotic evolution



Tahmasebi H, et al. Biomolecules 15(1): 93, 2025. PubMed: 39858487

Spread of Resistance

Healthcare-acquired infections (HAIs): infections that occur within 30 days of receiving health care

More than a half-million people get health care-associated infections (HAIs)


N Engl J Med 379(18): 1732-1744, 2018. PubMed: 30380384

Community-acquired infections (CAIs): infections that cannot be traced to healthcare

A recent study¹ of livestock manure identified:



- > 2000 known antibiotic-resistant genes (ARGs) resistant to 30 antibiotic classes
- ARG levels 2–18 X higher than in human feces, sewage, and soil.

¹Sci Adv 11(26): eadt8073, 2025. PubMed: 40577477




Spread of Resistance

Resistance can develop in multiple ways



Lancet Microbe 6(1): 100947, 2025. PubMed: 39305919

Antibacterial Discovery Pipeline



From Outterson, K. Health Affairs 2021

Meeting the Challenge: Global Efforts

- High-level political commitment
- Evidence-informed policy making
- One Health & multisectoral collaboration

WHO Three Global Priorities for Innovation:

- Severe multidrug-resistant (MDR) gram-negative infections
- Severe gram-positive infections in immunosuppressed and critically ill patients
- Bacterial meningitis, including penicillin- and cephalosporin-resistant community-acquired infections.



INNOVATIVE TECHNOLOGIES

- Novel chemistry
- Antimicrobial peptides
- CRISPR-Cas3
- AI/ML



ALTERNATIVES TO ANTIBIOTICS

- Phage therapy
- Immune modulating agents
- Monoclonal antibodies
- Microbiome



DIAGNOSTICS & SURVEILLANCE

- Rapid phenotypic testing

Meeting the Challenge: ATCC® Offerings



The most comprehensive, fully authenticated collection:


- 80,000+ bacteria, fungi, viruses, and protozoa
- 8,700+ microbial type strains
- 1,000+ derivatives, including genomic and synthetic nucleic acids
- 7,000+ microbial genomes
- 4,000+ cell lines and associated media and reagents

Explore our microbial products




Microbial Diversity

Bacteriology




1226 Genera

Mycology



1864 Genera

Virology




150 Genera

Brand recognition:

- Organizations and regulatory agencies specify ATCC® cultures - USP, ISO, FDA, CLSI, USDA, ASTM, AOAC, etc.
- Over 475 reference strains recommended for use in quality control

ATCC has Tools Supporting the Antibiotic Discovery Pipeline



Natural product producers and bacteriophages

Reporter cells

QC strains

AMR panels

Cell lines, hepatocytes, and organoids

Meeting the Challenge: ATCC® Offerings



GFP-tagged strains for antibiotic discovery



ATCC ID	Name	Reporter
25922GFP™	<i>Escherichia coli</i>	GFP
14028GFP™	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium	GFP
12022GFP™	<i>Shigella flexneri</i>	GFP
10145GFP™	<i>Pseudomonas aeruginosa</i>	GFP
15692GFP™	<i>Pseudomonas aeruginosa</i>	GFP

Case Study:

Use of 25922GFP™ to evaluate surface antibacterial activity



Meeting the Challenge: ATCC® Offerings

ATCC® strains are cited as QC controls in Clinical and Laboratory and Standards Institute (CLSI) AST methods



CLSI #	Method	ATCC® catalog number																										
		25922™	27853™	29212™	29213™	49619™	25923™	35218™	33186™	29741™	43055™	70057™	750™	6258™	22019™	24433™	90018™	90028™	MYA-3626™	MYA-3629™	MYA-3630™	MYA-3631™	MYA-3633™	MYA-3634™	MYA-4438™	MYA-3349™		
M02	Performance standards for disk Susceptibility tests																											
M07	Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically																											
M11	Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria																											
M23	Development of In Vitro Susceptibility Test Methods, Breakpoints, and Quality Control Parameters																											
M24	Performance Standards for Susceptibility Testing of <i>Mycobacteria</i> , <i>Nocardia</i> spp., and Other Aerobic Actinomycetes																											
M27	Performance Standards for Antifungal Susceptibility Testing of Yeasts																											
M38	Performance Standards for Antifungal Susceptibility Testing of Filamentous Fungi																											

Coming soon: Two new CLSI controls for Cefiderocol and Cefiderocol-xeruborbactam combination


Meeting the Challenge: ATCC® Offerings

Priority Antimicrobial-Resistant Strains



Selection

- Recent, clinically relevant isolates
- Prevalent pathogenic species
- Multidrug-resistant and extensively drug-resistant strains



Analysis

- *De novo* genome sequencing and annotation
- Identification of antimicrobial resistance genes
- Evaluation of antimicrobial resistance and susceptibility



Collection

- Strains authenticated and categorized by ATCC
- Genotypically, phenotypically, and functionally characterized

We are sourcing high-priority isolates to add to the collection

Multidrug-resistant
Pseudomonas aeruginosa from contaminated eye drops


BAA-3443™

BAA-3444™



Meeting the Challenge: ATCC® Offerings

WHO bacterial priority pathogens list, 2024



WHO Priority	Species	AMR	No of ATCC® strains
Critical	<i>Acinetobacter baumannii</i>	Carbapenem resistant	16
Critical	<i>Klebsiella pneumoniae</i>	Third-generation cephalosporin resistant	17
Critical	<i>Escherichia coli</i>	Third-generation cephalosporin resistant	19
Critical	<i>Proteus</i> spp.	Third-generation cephalosporin resistant	3
Critical	<i>Citrobacter</i> spp.	Third-generation cephalosporin resistant	8
Critical	<i>Serratia</i> spp.	Third-generation cephalosporin resistant	2
Critical	<i>Klebsiella pneumoniae</i>	Carbapenem resistant	6
Critical	<i>Escherichia coli</i>	Carbapenem resistant	12
Critical	<i>Enterobacter</i> spp.	Carbapenem resistant	7
Critical	<i>Citrobacter</i> spp.	Carbapenem resistant	2
Critical	<i>Mycobacterium tuberculosis</i>	Rifampicin resistant	2
High	<i>Enterococcus faecium</i>	Vancomycin resistant	9
High	<i>Pseudomonas aeruginosa</i>	Carbapenem resistant	15
High	<i>Neisseria gonorrhoeae</i>	Fluoroquinolone resistant	5
High	<i>Staphylococcus aureus</i>	Methicillin resistant	69
Medium	Group A streptococci	Macrolide resistant	5
Medium	<i>Streptococcus pneumoniae</i>	Macrolide resistant	28
Medium	<i>Haemophilus influenzae</i>	Ampicillin resistant	2

Meeting the Challenge: ATCC® Offerings

Priority antimicrobial-resistant strains



Acinetobacter baumannii Bouvet and Grimont

BAA-3252™

Download Genome

LEARN ABOUT THE ATCC GENOME PORTAL >

This drug-resistant bacterium was isolated in 2014 from the abscess of a 57-year-old male in Spain. This product can be used for all stages of the discovery and development process for novel antimicrobials and therapeutics, molecular-based detection assays, and updated sterility protocols.

Product category	Bacteria
Product type	Drug-resistant bacterium
Strain designation	1074318
Type strain	No
Genome sequenced strain	Yes
Isolation source	Abscess
Geographical isolation	Spain
Applications	Antimicrobial resistance research Bioinformatics Drug development
Product format	Frozen
Storage conditions	-80°C or colder

Buy Now

Price: **\$297.00 ea**

Discounts may be available for our fellow nonprofit organizations. [Login](#) to see your price.

Generally ships within 1-3 business days

Quantity

[Add to cart](#)

[Add to list](#)


Documentation

- Product sheet
- Certificate of analysis
- Safety data sheet
- Antimicrobial resistance report
- Characterization data**



Priority AMR Strains Brochure


www.atcc.org/AMR-brochure



Meeting the Challenge: ATCC® Offerings

Acinetobacter baumannii: Priority Antimicrobial-Resistant Strains




Antibiotics	Class	Antibiotic target	BAA-3278™	BAA-3301™	BAA-3275™	BAA-3320™	BAA-3257™	BAA-3276™	BAA-3300™	BAA-3302™	BAA-3252™	BAA-3282™	BAA-3283™	BAA-3338™	BAA-3311™	
Gentamicin	Aminoglycosides	30S	≥16 (R)	≥16 (R)	≤1 (S)	≤1 (S)	≥16 (R)	4 (S)	≥16 (R)	8 (I)	≤1 (S)	4 (S)	≤1 (S)	≥16 (R)	≤1 (S)	
Tobramycin			≥16 (R)	≥16 (R)	≤1 (S)	≤1 (S)	≥16 (R)	2 (S)	≥16 (R)	8 (I)	≤1 (S)	8 (I)	2 (S)	8 (I)	≤1 (S)	
Imipenem	Carbapenems	Cell wall	8 (I)	8 (I)	8 (I)	≥16 (R)	8 (I)	≥16 (R)	8 (I)	≥16 (R)	≥16 (R)	≥16 (R)	≥16 (R)	≥16 (R)	0.25 (S)	
Meropenem			≥16 (R)	4 (S)	≥16 (R)	≥16 (R)	4 (S)	≥16 (R)	≥16 (R)	8 (I)	≥16 (R)	≥16 (R)	≥16 (R)	≥16 (R)	≥16 (R)	0.25 (S)
Cefazolin	Third Generation Cephalosporin	Cell wall	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	
Cefepime			≥64 (R)	16 (I)	32 (R)	≥64 (R)	8 (S)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	32 (R)	≥64 (R)	2 (S)
Cefotaxime			≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	16 (I)
Ceftazidime			≥64 (R)	≥64 (R)	≥64 (R)	16 (I)	16 (I)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	4 (S)
Ceftriaxone			≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	≥64 (R)	16 (I)
TMP-sulfamethoxazole	Diaminopyrimidines	Folic acid metabolism	160 (R)	160 (R)	≥320 (R)	≤20 (S)	≥320 (R)	160 (R)	≤20 (S)	≥320 (R)	≥320 (R)	≥320 (R)	≥320 (R)	≥320 (R)	≤20 (S)	
Ciprofloxacin	Fluoroquinolones	DNA gyrase	≥4 (R)	≥4 (R)	2 (I)	1 (S)	≥4 (R)	≥4 (R)	≥4 (R)	≥4 (R)	≥4 (R)	≥4 (R)	≥4 (R)	≥4 (R)	0.25 (S)	
Levofloxacin			≥8 (R)	≥8 (R)	1 (S)	0.5 (S)	≥8 (R)	≥8 (R)	≥8 (R)	≥8 (R)	≥8 (R)	≥8 (R)	≥8 (R)	≥8 (R)	≥8 (R)	0.12 (S)
Ampicillin/Sultbactam	Penicillin/BLI	Cell wall	16 (I)	4 (S)	8 (S)	16 (I)	16 (I)	8 (S)	16 (I)	≥32 (R)	8 (S)	≤2 (S)	16 (I)	≥32 (R)	≤2 (S)	
Piperacillin/Tazaobactam			≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	≥128 (R)	8 (I)
Tetracycline	Tetracyclines	30S	2 (S)	2 (S)	≥16 (R)	2 (S)	2 (S)	4 (S)	2 (S)	2 (S)	≥16 (R)	≥16 (R)	≥16 (R)	4 (S)	≤1 (S)	
	Metadata	Isolation country	Honduras	Italy	France	Venezuela	Spain	USA	Singapore	Italy	Spain	USA	Taiwan	USA	Korea	
		Year of isolation	2006	2006	2006	2006	2004	2014	2004	2015	2014	2004	2006	2004	2006	
		Source of isolation	Catheters	Blood	Blood	Blood	Blood	Sputum	Trachea	Blood	Abscess	Wound	Sputum	Blood	Blood	
		Patient Gender	Female	Female	Male	Male	Female	Male	Male	Male	Male	Female	Female	Female	Female	Female
		Patient Age (Yr)	<1	53	84	19	68	75	71	41	57	51	80	764	73	

Priority AMR Strains Brochure
www.atcc.org/AMR-brochure

Meeting the Challenge: ATCC® Offerings

New bacterial isolates coming soon




WHO Priority	Species	AMR	No of ATCC® strains	Why
Critical	<i>Enterobacter</i> spp.	Third-generation cephalosporin resistant	8	Recent isolates
Critical	<i>Escherichia coli</i>	Third-generation cephalosporin resistant	3	Recent isolates
Critical	<i>Morganella morganii</i>	Third-generation cephalosporin resistant	3	Gap
Critical	<i>Pluralibacter</i> spp.	Third-generation cephalosporin resistant	1	Gap
Critical	<i>Enterobacter</i> spp.	Carbapenem resistant	3	Recent isolates
Critical	<i>Morganella morganii</i>	Carbapenem resistant	2	Gap
Critical	<i>Klebsiella pneumoniae</i>	Carbapenem resistant	1	Recent isolates
Critical	<i>Shigella</i> spp.	Fluoroquinolone resistant	1	Gap

- *Mycobacterium* strains resistant to known treatments developed in an attenuated mycobacterial strain (TB vaccine strain Bacillus Calmette-Guérin, or BCG) containing defined drug resistance mutations.
- New two new CLSI reference materials for Cefiderocol and Cefiderocol-xeruborbactam

Meeting the Challenge: ATCC® Offerings

WHO fungal priority pathogens 2022 list



WHO Priority	Species	No of ATCC® strains
Critical	<i>Cryptococcus neoformans</i>	134
Critical	<i>Candida auris</i>	8
Critical	<i>Aspergillus fumigatus</i>	101
Critical	<i>Candida albicans</i>	227
High	<i>Nakaseomyces glabrata</i>	30
High	<i>Histoplasma</i> spp.	62
High	Eumycetoma causing agents	59
High	Mucorales	287
High	<i>Fusarium</i> spp.	1125
High	<i>Candida tropicalis</i>	67
High	<i>Candida parapsilosis</i>	41
Medium	<i>Paracoccidioides</i> spp.	31
Medium	<i>Cryptococcus gattii</i>	26
Medium	<i>Pichia kudriavzevii</i>	20
Medium	<i>Talaromyces marneffei</i>	19
Medium	<i>Scedosporium</i> spp.	18
Medium	<i>Coccidioides</i> spp.	11
Medium	<i>Lomentospora prolificans</i>	7
Medium	<i>Pneumocystis jirovecii</i>	1

Meeting the Challenge: ATCC® Offerings

New fungal pathogens coming soon



Strain	# of Isolates	Amphotericin B	Anidulafungin	Caspofungin	Fluconazole	Micafungin	Itraconazole	Terbinafine	Gap in Collection
<i>Candida auris</i>	5	Red	Green	Red	Red	Green			Gap, fill caspofungin and amphotericin resistance
<i>Candida auris</i>	15	Red	Green	Green	Red	Green			Gap, fill caspofungin resistance
<i>Candida auris</i>	11	Green	Green	Green	Red	Green			New isolates
<i>Candida auris</i>	1	Green	Red	Red	Red	Red			Gap, fill Anidulafungin resistance
<i>T. indotineae</i> ¹	5						Green	Red	Gap, new species
<i>T. indotineae</i>	2						Red	Green	Gap, new species
<i>T. indotineae</i>	2						Red	Red	Gap, new species
<i>T. indotineae</i>	3						Green	Green	Gap, new species

¹*Trichophyton indotineae*

ATCC® Genome Portal

The only authenticated reference genome database for ATCC microbes

The ATCC® Genome Portal (AGP) is a rapidly growing ISO 9001-compliant database of high-quality reference genomes from authenticated microbial strains in the ATCC collection. Customers can easily access and download meticulously curated whole-genome assemblies for purchased strains and *Supporting Members* have full access to the AGP.




7,000+

Available microbial reference genomes
as of June 2026

- Download ATCC® microbial genomes.
- View ClinVar variants for ATCC® cell-lines.
- Search by sequence and annotations.
- Stratify taxonomy, isolation source, catalog number, type strain status, strain name, biosafety level, etc.
- View genome assembly statistics and quality metrics.
- Find relatedness of ATCC® genomes.

New genomes are released at the end of every quarter.



Learn more about the
ATCC® Genome Portal

Vial to File


Physical Repository

- Strains
- Derivatives
- Standards
- Reference materials




Authenticated Reference Data

- Sequencing data
- Assembled genomes
- Annotated genes



Drive scientific advancement by providing the scientific community with high-quality, annotated whole-genome sequence (WGS) information to complement ATCC's biological materials.

Authenticated Physical Material Coupled with Reference-Quality Genome Sequences



- Fully traceable and authenticated to ATCC® materials.
- All genome assemblies produced in-house at ATCC® in an ISO certified laboratory.

ATCC® Genome Portal – Ways to visit the AGP



Home > Microbe Products > Bacteriology and Archaea > BAA-3330

Streptococcus pneumoniae (Klein) Chester

BAA-3330™

[Download Genome](#)

[LEARN ABOUT THE ATCC GENOME PORTAL >](#)

This drug-resistant bacterium was isolated in 2013 from the respiratory bronchoalveolar lavage of a 22-year-old male in Italy. This product can be used for all stages of the discovery and development process for novel antimicrobials and therapeutics, molecular-based detection assays, and updated sterility protocols.

Product category	Bacteria
Product type	Drug-resistant bacterium
Strain designation	978059
Type strain	No
Genome sequenced strain	Yes

Buy Now

Price: \$304.00 ea

Discounts may be available for our fellow nonprofit organizations. [Login](#) to see your price.

Generally ships within 1-3 business days

Quantity

[Add to cart](#)

[Add to list](#)

Documentation

- [Product sheet](#)
- [Certificate of analysis](#)
- [Safety data sheet](#)
- [Antimicrobial resistance report](#)
- [Characterization data](#)

Careers Support United States

Log In Create a Profile Order Status Quick Order

Home > Applications >

Discover

7,000 ge
3,000 tra

[Visit the ATCC Genome](#)

Microbe Products

- By Application
- [ATCC Genome Portal](#)
- Virology
- Bacteriology and Archaea
- Mycology
- Type Strains
- Collections and Projects
- Protistology

MicroQuant™

Streamline microbial QC testing with precisely quantitated, ready-to-use controls

[LEARN MORE >](#)

Authentic, standardized omics data for ATCC microbes and cell lines

Researchers need rigorously validated, traceable genomic data to make reliable scientific correlations—but publicly available sequences often lack authenticity, quality, and provenance. To close this gap, ATCC developed standardized sequencing, assembly, and annotation pipelines that deliver high-quality reference genomes from our authenticated microbial strains, along with whole-exome and RNA-seq data from our cell lines. These datasets are accessible through the [ATCC Genome Portal](#)—a rapidly growing, ISO 9001-compliant, cloud-based platform offering authenticated NGS data, comprehensive metadata, and advanced tools for reproducible analyses.

[Request a demo](#) [Visit the ATCC Genome Portal](#) [Explore our resources](#) [Watch our tutorial](#)



ATCC® Genome Portal
www.atcc.org/genomeportal

ATCC® Genome Portal – Browse for Strains / Data



- HOME
- MICROBE PRODUCTS
- CELL PRODUCTS
- DOCUMENTATION

JPETRONE@ATCC.ORG

Welcome to the ATCC Genome Portal

The only authenticated reference omics database for ATCC bi...

[VIEW ALL MICROBES >](#)

[VIEW ALL CELL LINES >](#)

ATCC Genome Portal Qtr 2 Updates - June 03, 2025

- Released new datasets for **49** human and **2** organoid cell lines
- Added genomes for **225** Bacteria, **7** Viruses, **16** Fungi, and **2** Protists
- Published **32** genomes supporting the FDA's **GenomeTrakr** initiative
 - To query, enter "**GenomeTrakr**" in the **advanced search feature** and select **GenomeTrakr** in the **Tag** filter.
- Expanded to include **4.5K** epigenome and **5.3K** bacterial resistance datasets

- HOME
- MICROBE PRODUCTS
- CELL PRODUCTS
- DOCUMENTATION

Genomes

[All Microbes](#)

0 out of 7000 selected (Select All / Deselect All)

<input type="checkbox"/>	Taxonomic Name	ATCC Product Page	Strain Alias	Tags	Date Published	Length
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 49176™	SC10 [CIP 103242]	Methylation Data Type Strain	December 12, 2022	2.0 Mb
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 700209™	CIP 103893 [PE7]	Methylation Data	April 29, 2024	2.0 Mb
<input type="checkbox"/>	<i>Acetivibrio aldrichii</i>	ATCC® 49358™	P-1 [OGI 112]	Methylation Data Type Strain	September 25, 2024	6.4 Mb
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 33288™	CD2 [NRCC 2248]	Methylation Data Type Strain	November 26, 2024	6.3 Mb
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 35928™	BAS [NRCC 2936]	Methylation Data Type Strain	March 27, 2025	6.3 Mb
<input type="checkbox"/>	<i>Acetivibrio ethanoligignens</i>	ATCC® 33324™	77-6	Methylation Data Type Strain	June 3, 2024	4.1 Mb
<input type="checkbox"/>	<i>Acetivibrio thermocellus</i>	ATCC® 27405™	VPI 7372 [157]	Type Strain	August 27, 2019	3.8 Mb
<input type="checkbox"/>	<i>Acetoanaerobium noterae</i>	ATCC® 35199™	NOT-3	Methylation Data Type Strain	June 11, 2025	2.9 Mb
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 15973™	NCIB 8621 [Delft L. 40, J. Frateur 4 (Collection I)]	Methylation Data Type Strain	September 29, 2020	3.7 Mb
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 23746™	NCIB 8554	Methylation Data	January 28, 2021	3.7 Mb
<input type="checkbox"/>	<i>Acetobacterium sp</i>	ATCC® BAA-990™	SyrA5	Methylation Data	March 5, 2024	4.1 Mb
<input type="checkbox"/>	<i>Acetobacter orleanensis</i>	ATCC® 12876™	NCIB 8622 [DSM 4492, IFO 13752, JCM 7639, L 1027, LMG 1583, NCCB 31003]	Methylation Data Type Strain	June 11, 2025	3.1 Mb

ATCC® Genome Portal – Reference Genome Pages

Salmonella enterica subsp. *enterica* serovar Bareilly (ATCC® BAA-3403™)



Overview page

ATCC | HOME | MICROBE PRODUCTS | CELL PRODUCTS | DOCUMENTATION | Quick Search... | JPETRONE@ATCC.ORG

Salmonella enterica subsp. enterica > **Salmonella enterica subsp. enterica serovar Bareilly**

Salmonella enterica subsp. enterica serovar Bareilly (ATCC® BAA-3403™)

Overview | Genome Browser | Related Genomes | Quality Control

DOWNLOAD DATA | RUN DISCREPANCY REPORT

Organism Summary		Assembly Summary	
Catalog Number	ATCC® BAA-3403™	Date Published	September 16, 2025
Tags	GenomeTrakr Methylation Data	Length	4,808,140 nt
Type Strain	No	Sequencing Technology	Illumina + Oxford Nanopore Hybrid Assembly
		Number of Contigs	2 (All Circularized)
		N50	4,726,325 nt
		%GC	52.21%

Annotations Summary	
Number of CDS	4534
Number of Hypothetical Proteins	184
Number of tRNA	82
Number of 5s rRNA	8
Number of 16s rRNA	7
Number of 23s rRNA	7

Genome Annotations Browser

ATCC | HOME | MICROBE PRODUCTS | CELL PRODUCTS | DOCUMENTATION | Quick Search... | JPETRONE@ATCC.ORG

Salmonella enterica subsp. enterica > **Salmonella enterica subsp. enterica serovar Bareilly**

Salmonella enterica subsp. enterica serovar Bareilly (ATCC® BAA-3403™)

Overview | **Genome Browser** | Related Genomes | Quality Control

Annotation Legend

Include Hypothetical Proteins | Display All Genes | Filter annotations | Download Table CSV

Contig	Start	End	Name	Protein Product	EC Number	Type	RefSeq ID	Jump
1	1	1401	<i>dnaA</i>	chromosomal replication initiator protein DnaA		CDS	NP_462738.1	Jump
1	1406	2506	<i>dnaN</i>	DNA polymerase III subunit beta	2.7.7.7	CDS	WP_006177590.1	Jump
1	2654	3727	<i>recF</i>	DNA replication/repair protein RecF		CDS	WP_005121479.1	Jump
1	3756	6170	<i>gyrB</i>	DNA topoisomerase (ATP-hydrolyzing) subunit B	5.6.2.2	CDS	WP_005121480.1	Jump
1	7085	6189		LysR family transcriptional regulator		CDS	NP_462734.1	Jump
1	8393	7200		mandelate racemase/muconate lactonizing enzyme family protein		CDS	NP_462733.1	Jump
1	9669	8407		MFS transporter		CDS	NP_462732.1	Jump
1	9974	10819	<i>yidA</i>	sugar-phosphatase	3.1.3.23	CDS	NP_462731.1	Jump
1	11080	11769	<i>dgoR</i>	D-galactonate utilization transcriptional regulator DgoR		CDS	NP_462730.1	Jump
1	11766	12644		2-dehydro-3-deoxygalactonokinase		CDS	NP_462729.1	Jump

ATCC® Genome Portal – Reference Genome Pages

Salmonella enterica subsp. *enterica* serovar Bareilly (ATCC® BAA-3403™)



Related Genomes

ATCC® HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search... J.PETRONI@ATCC.ORG

Salmonella enterica subsp. enterica > Salmonella enterica subsp. enterica serovar Bareilly

Salmonella enterica subsp. enterica serovar Bareilly (ATCC® BAA-3403™)

Overview Genome Browser Related Genomes Quality Control

Most similar genomes

The following genomes have the greatest genomic similarity to this one (≥95% average nucleotide identity):

All Microbes My Favorite Microbes

DOWNLOAD DATA 0 out of 81 selected (Select All / Deselect All) Type for Advanced Search...

<input type="checkbox"/>	Taxonomic Name	ATCC Product Page	Similarity	Strain Alias	Tags	Date Published	Length
<input type="checkbox"/>	★ <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Bareilly strain GTVSS-001	ATCC® BAA-3401™	96.42%		GenomeTrakr	November 21, 2025	4.8 Mb
<input type="checkbox"/>	★ <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Bareilly	ATCC® 9115™	96.33%		Methylation Data	October 29, 2021	4.8 Mb
<input type="checkbox"/>	★ <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Saintpaul	ATCC® 9712™	95.46%	127	Methylation Data	February 24, 2023	4.7 Mb
<input type="checkbox"/>	★ <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Thompson	ATCC® BAA-3141™	95.46%			August 1, 2022	4.7 Mb
<input type="checkbox"/>	★ <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Stanley	ATCC® 7308™	95.46%		Methylation Data	October 28, 2020	4.7 Mb
<input type="checkbox"/>	★ <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Thompson	ATCC® 8391™	95.46%	2988		December 15, 2021	4.7 Mb

Quality Control

ATCC® HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search... J.PETRONI@ATCC.ORG

Salmonella enterica subsp. enterica > Salmonella enterica subsp. enterica serovar Bareilly

Salmonella enterica subsp. enterica serovar Bareilly (ATCC® BAA-3403™)

Overview Genome Browser Related Genomes Quality Control

Sequencing Quality Control

Sequencing Quality Control Metrics	
Number of trimmed reads	5,285,738
Median Q score, all bases	40
Percent of median Q scores per cycle greater than 25	1%
Ambiguous content (% N bases)	0%

Assembly Quality Control

Assembly Quality Control Metrics	
Estimated genome completeness	99.64%
Estimated genome contamination	0.81%
Illumina depth of coverage	163.856x
ONT depth of coverage	61.48x


Notes

This bacterium was assembled, annotated, and curated using Oatmeal v3.0. See our technical documentation for further details: <https://docs.onecodex.com/en/articles/5357065-assembly-qc-and-annotation-pipeline-versions>. This item is part of the GenomeTrakr Validation Strain Set (GTVSS) as foodborne and environmental pathogens controls to support validation and performance testing of whole-genome sequencing (WGS) workflows. Please see the MP-44 - GTVSS product page for more information: <https://www.atcc.org/products/mp-44>

ATCC[®] Genome Portal – Bacterial Resistome Pipeline



- Hybrid assembly and ILM polishing gives confidence to SNP-based AMR.
- Resistome datasets are uploaded with every new bacterial genome deposition.
- Resistome packages also updated yearly with PGAP and AMR tool updates.
- VITEK 2 and Antibigram data also made available from AGP and JSON format.



ATCC® Genome Portal – Ways to view AMR

Via Product page / GUI



ATCC®

HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search... Become a Supporting Member LOG IN

Genomes

All Microbes

DOWNLOAD DATA 0 out of 6750 selected (Select All / Deselect All).

<input type="checkbox"/>	Taxonomic Name ^	ATCC Product Page	Strain Alias	Tags
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 49176™	SC10 [CIP 103242]	Methylation Data Type
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 700209™	CIP 103893 [PE7]	Methylation Data
<input type="checkbox"/>	<i>Acetivibrio aldrichii</i>	ATCC® 49358™	P-1 [OGI 112]	Methylation Data Type
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 33288™	CD2 [NRCC 2248]	Methylation Data Type
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 35928™	BAS [NRCC 2936]	Type Strain
<input type="checkbox"/>	<i>Acetivibrio ethanolignens</i>	ATCC® 33324™	77-6	Methylation Data Type
<input type="checkbox"/>	<i>Acetivibrio thermocellus</i>	ATCC® 27405™	VPI 7372 [157]	Type Strain
<input type="checkbox"/>	<i>Acetoanaerobium noterae</i>	ATCC® 35199™	NOT-3	Methylation Data Type
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 15973™	NCIB 8621 [Delft L. 40, J. Frateur 4 (Collection I)]	Methylation Data Type
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 23746™	NCIB 8554	Methylation Data
<input type="checkbox"/>	<i>Acetobacterium sp</i>	ATCC® BAA-990™	SyrA5	Methylation Data

BAA-3330

Name
BAA-3330

Catalog Number
BAA-3330

Strain Alias
BAA-3330

Tag
BAA-3330

Taxonomy
BAA-3330

Isolation
BAA-3330

Type Strain
BAA-3330

Biosafety Level
BAA-3330

Gene
BAA-3330

March 9, 2024 4.1 Mb

ATCC® Genome Portal – Ways to view AMR

Via Product page / GUI



The screenshot shows the ATCC Genome Portal interface. At the top, there is a navigation bar with links for HOME, MICROBE PRODUCTS, CELL PRODUCTS, and DOCUMENTATION, along with a search bar and a 'Become a Supporting Member' button. The main content area is titled 'Genomes' and shows a list of products. The first product is *Streptococcus pneumoniae*, which is highlighted. The 'Tags' column for this product is circled in blue, showing 'AMR' and 'Methylation Data'. A blue arrow points from the 'New' button to the search bar. The interface also includes a 'DOWNLOAD DATA' button and a 'Displaying 1 Product' message.

<input type="checkbox"/>	Taxonomic Name ▲	ATCC Product Page ⇅	Strain Alias	Tags	Date Published ⇅	Length
<input type="checkbox"/>	<i>Streptococcus pneumoniae</i>	ATCC® BAA-3330™ ↗	978059	AMR Methylation Data	June 11, 2025	2.1 Mb

ATCC® Genome Portal – Ways to view AMR

Via Product page / GUI



ATCC

HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search... *i*

Become a Supporting Member *LOG IN*

Streptococcus > Streptococcus pneumoniae

Streptococcus pneumoniae (ATCC® BAA-3330™) *↗*

Annotation Legend

Overview Genome Browser Related Genomes Quality Control

pbp2X

Include Hypothetical Proteins Display AMR Genes AMR Filters

Download Table CSV

Contig	Start	End	Name	Protein Product	EC Number	Type	RefSeq ID	Jump
1	1072482	1072054	<i>aac(6')</i>	aminoglycoside 6'-N-acetyltransferase	2.3.1.82	CDS	WP_013903432.1	<i>↗</i>
1	1153448	1152249		multidrug efflux MFS transporter		CDS	WP_000136121.1	<i>↗</i>
1	1511931	1509874	<i>pbp2b</i>	penicillin-binding protein PBP2B		CDS	WP_001829432.1	<i>↗</i>
1	1764181	1764918	<i>erm(B)</i>	23S rRNA (adenine(2058)-N(6))-methyltransferase Erm(B)		CDS	WP_011100845.1	<i>↗</i>
1	1774833	1776752	<i>tet(M)</i>	tetracycline resistance ribosomal protection protein Tet(M)		CDS	WP_011117014.1	<i>↗</i>
1	1779941	1778478	<i>msr(D)</i>	ABC-F type ribosomal protection protein Msr(D)		CDS	WP_000420317.1	<i>↗</i>
1	1781278	1780061	<i>mef(A)</i>	macrolide efflux MFS transporter Mef(A)		CDS	WP_000417520.1	<i>↗</i>
1	1953479	1951713	<i>patB</i>	multidrug efflux ABC transporter subunit PatB		CDS	WP_000859866.1	<i>↗</i>
1	1955965	1954271	<i>patA</i>	multidrug efflux ABC transporter subunit PatA		CDS	WP_000908145.1	<i>↗</i>

ATCC[®] Genome Portal – AMR Gene Search: Gene Name and/or Sequence



Sequence Searching

ATCC® Genome Portal – AMR Gene Search: Gene Name and/or Sequence



ATCC®

HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search... Become a Supporting Member LOG IN

Genomes

All Microbes

DOWNLOAD DATA 0 out of 7000 selected (Select All / Deselect All)

<input type="checkbox"/>	Taxonomic Name ▲	ATCC Product Page	Strain Alias	Tag	Length
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 49176™	SC10 [CIP 103242]	Methylation Data	2.0 Mb
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 700209™	CIP 103893 [PE7]	Methylation Data	2.0 Mb
<input type="checkbox"/>	<i>Acetivibrio aldrichii</i>	ATCC® 49358™	P-1 [OGI 112]	Methylation Data	6.4 Mb
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 33288™	CD2 [NRCC 2248]	Methylation Data	6.3 Mb
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 35928™	BAS [NRCC 2936]	Methylation Data Type Strain	March 27, 2025 6.3 Mb
<input type="checkbox"/>	<i>Acetivibrio ethanolgignens</i>	ATCC® 33324™	77-6	Methylation Data Type Strain	June 3, 2024 4.1 Mb
<input type="checkbox"/>	<i>Acetivibrio thermocellus</i>	ATCC® 27405™	VPI 7372 [157]	Type Strain	August 27, 2019 3.8 Mb
<input type="checkbox"/>	<i>Acetoanaerobium noterae</i>	ATCC® 35199™	NOT-3	Methylation Data Type Strain	June 11, 2025 2.9 Mb
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 15973™	NCIB 8621 [Delft L. 40, J. Frateur 4 (Collection I)]	Methylation Data Type Strain	September 29, 2020 3.7 Mb
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 23746™	NCIB 8554	Methylation Data	January 28, 2021 3.7 Mb
<input type="checkbox"/>	<i>Acetobacterium sp.</i>	ATCC® BAA-990™	SvrA5	Methylation Data	March 5, 2024 4.1 Mb


ATTTATTGCTATTTA CCGCGGCTTTTATTGAGCTTGAAAGATAAATAA

Search for genomes containing the sequence (minimum 31 bases):

ATGGAAAATAAAAAATCATCAACAAGAAAATTTAAGAGTAC2TATC
 AATCACTGGTTAACTCAGCACGAATATTGTTTGTGAAAAAGGCTA
 TCAAGCTGTTCAATAGATGAGATCTCGGAAAAGCGTTGGTGACC
 AAAGGTGCCTTTTATCATCCTTTAAAAATAAAAAACAATTACTCAG
 TGCCTGTATAAGCAGCAATTAATTATGATTGATGCCTACATCACAA
 CAAAACTGATTTAACAATGGTTGGTCTGCCTTAGAAAGTATATT
 TGAACATATCTTGATTATATTATTGATAATAATAAAAAACCTTATCC
 TATCCAAGAAGTGATGCCTATCATTGGTTGGAATGAACCTGAAAAA
 ATTAGCCTTGAATACATTACTGGTAAGGTAACGCCATTGTCAGCA
 AATTGATCCAAGAGAACCAACTAAAGCTTATGATGATGATGTGCT
 TAAAACTTACTCAATGGCTGGTTTATGCATATCGCAATACATGCGA
 AAAACCTAAAAGAGCTTGCCGATAAAAAAGCCAATTTATTGCTAT
 TTACCGCGGCTTTTATTGAGCTTGAAAGATAAATAA

ATCC® Genome Portal – AMR Gene Search: Gene Name and/or Sequence





[HOME](#)
[MICROBE PRODUCTS](#)
[CELL PRODUCTS](#)
[DOCUMENTATION](#)

[Become a Supporting Member](#)

[LOG IN](#)

Genomes

[All Microbes](#)

Nucleotide search for 594bp query i

```

ATGGAAAATAAAAATCATCAACAAGAAAATTTAAGAGTACCTATCAATCACTGGTTAACTCAGCACGAATATTGTTTGTGAAAAAGGCTATCAAGCTGTTTCAATAGATGAGATCTCGGGAAAAGCGTTGGTGACCAAAGGTGCCTTTTATCATCACTT
TAAAAATAAAAAACAATTAAGAGTACCTATCAATCACTGGTTAACTCAGCACGAATATTGTTTGTGAAAAAGGCTATCAAGCTGTTTCAATAGATGAGATCTCGGGAAAAGCGTTGGTGACCAAAGGTGCCTTTTATCATCACTT
TCCCTATCCAAGAAGTGATGCCATCATTGGTTGGAATGAACCTGAAAAAATTAGCCTTGAATACATTACTGGTAAGGTAACGCCATTGTGAGCAAAATGATCCAAGAGAACCAACTTAAAGCTTATGATGATGATGTGCTTAAAAACTTACTCAATGGC
TGGTTTATGCATATCGCAATACATGCGAAAAACCTAAAAGAGCTTCCGCGATAAAAAAGGCCAATTTATTGCTATTTACCGCGGCTTTTATTGAGCTTGAAGATAAATAA
                    
```

↓ [DOWNLOAD DATA](#)
0 out of 26 selected (Select All / Deselect All).

Type for Advanced Search...

☰
🏠

<input type="checkbox"/>	Taxonomic Name ↕	ATCC Product Page ↕	Similarity ▼	Strain Alias	Tags	Date Published ↕	Length
<input type="checkbox"/>	<i>Shigella dysenteriae</i>	ATCC® 49552™ ↗	594 bases matched (100%)	CDC 3010-85	Methylation Data	March 24, 2025	5.0 Mb
<input type="checkbox"/>	<i>Klebsiella pneumoniae subsp pneumoniae</i>	ATCC® 27799™ ↗	594 bases matched (100%)	1296	Methylation Data	April 28, 2023	6.0 Mb
<input type="checkbox"/>	<i>Escherichia coli</i>	ATCC® BAA-2778™ ↗	594 bases matched (100%)	1123005	AMR	February 26, 2020	5.0 Mb
<input type="checkbox"/>	<i>Moellerella wisconsinensis</i>	ATCC® 35620™ ↗	594 bases matched (100%)	CDC 2897-78	Methylation Data	December 12, 2025	3.5 Mb
<input type="checkbox"/>	<i>Escherichia coli</i>	ATCC® 27065™ ↗	594 bases matched (100%)	2395 [NCIB 12415]	Methylation Data	September 16, 2025	4.7 Mb
<input type="checkbox"/>	<i>Escherichia coli</i>	ATCC® BAA-1161™ ↗	594 bases matched (100%)	UMN 026	Methylation Data	April 2, 2021	5.4 Mb

ATCC® Genome Portal – AMR Gene Search: Gene Name and/or Sequence



Annotation Searching

ATCC® Genome Portal – AMR Gene Search: Gene Name and/or Sequence



ATCC®

HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search... Become a Supporting Member LOG IN

Genomes

All Microbes

DOWNLOAD DATA 0 out of 7000 selected (Select All / Deselect All)

<input type="checkbox"/>	Taxonomic Name ▲	ATCC Product Page ↗	Strain Alias	Tag	Length
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 49176™ ↗	SC10 [CIP 103242]	Me	2.0 Mb
<input type="checkbox"/>	<i>Abiotrophia defectiva</i>	ATCC® 700209™ ↗	CIP 103893 [PE7]	Me	2.0 Mb
<input type="checkbox"/>	<i>Acetivibrio aldrichii</i>	ATCC® 49358™ ↗	P-1 [OGI 112]	Me	6.4 Mb
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 33288™ ↗	CD2 [NRCC 2248]	Me	6.3 Mb
<input type="checkbox"/>	<i>Acetivibrio cellulolyticus</i>	ATCC® 35928™ ↗	BAS [NRCC 2936]	Me	6.3 Mb
<input type="checkbox"/>	<i>Acetivibrio ethanolignens</i>	ATCC® 33324™ ↗	77-6	Me	4.1 Mb
<input type="checkbox"/>	<i>Acetivibrio thermocellus</i>	ATCC® 27405™ ↗	VPI 7372 [157]	Typ	3.8 Mb
<input type="checkbox"/>	<i>Acetoanaerobium noterae</i>	ATCC® 35199™ ↗	NOT-3	Me	2.9 Mb
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 15973™ ↗	NCIB 8621 [Delft L. 40, J. Frateur 4 (Collection I)]	Me	3.7 Mb
<input type="checkbox"/>	<i>Acetobacter aceti</i>	ATCC® 23746™ ↗	NCIB 8554	Me	3.7 Mb
<input type="checkbox"/>	<i>Acetobacterium sp</i>	ATCC® BAA-990™ ↗	SvrA5	Methylation Data	March 5, 2024 4.1 Mb

blaM

Name

blaM

Catalog Number

blaM

Strain Alias

blaM

Tag

blaM

Taxonomy

blaM

Isolation

blaM

Type Strain

blaM

Biosafety Level

blaM

Gene

blaM

ATCC® Genome Portal – AMR Gene Search: Gene Name and/or Sequence



ATCC®

HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search... Become a Supporting Member LOG IN

... > *Aeromonas* > ***Aeromonas caviae***

Aeromonas caviae (ATCC® 15468™)

Overview Genome Browser Related Genomes Quality Control Annotation Legend

Include Hypothetical Proteins Display AMR Genes AMR Filters Q blaM Download Table CSV

Contig	Start	End	Name	Protein Product	EC Number	Type	RefSeq ID	Jump
1	3340859	3339708	<i>blaMOX</i>	CMY-1/MOX family class C beta-lactamase	3.5.2.6	CDS	WP_043155783.1	▶


Powered by ONE CODEX

Discover the ATCC Genome Portal Data Use Agreement Need Help? Contact us at ATCC Support

ATCC® Genome Portal – Ways to view AMR

View / Download Phenotypical Data



HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION ? Become a Supporting Member LOG IN

... > *Streptococcus* > *Streptococcus pneumoniae*

Streptococcus pneumoniae (ATCC® BAA-3330™) [↗](#)

Overview Genome Browser Related Genomes **Quality Control**

Sequencing Quality Control

Sequencing Quality Control Metrics	
Number of trimmed reads	1,908,100
Median Q score, all bases	40
Percent of median Q scores per cycle greater than 25	1%
Ambiguous content (% N bases)	0%

Assembly Quality Control

Assembly Quality Control Metrics ?	
Estimated genome completeness	99.62%
Estimated genome contamination	0.31%
Illumina depth of coverage	238.621x
ONT depth of coverage	864.051x

Notes

This bacterium was assembled, annotated, and curated using Oatmeal v3.0. See our technical documentation for further details: <https://docs.onecodex.com/en/articles/5357065-assembly-qc-and-annotation-pipeline-versions>. This item includes additional phenotypic antimicrobial resistance (AMR) data. You may download the Vitek2 antibiotic susceptibility report here: <https://www.atcc.org/-/media/product-assets/documents/antimicrobial-susceptibility-reports/BAA-3330-antibiotic-susceptibility-profile.pdf>

ATCC® Genome Portal – Ways to view AMR

View / Download Phenotypical Data



Antibiotic Susceptibility of *Streptococcus pneumoniae* Strain 978059 ATCC® BAA-3330™

Antimicrobial	MIC ^a	Interpretation ^b	Antimicrobial	MIC	Interpretation
Benzylpenicillin (Meningitis)	≥8	R	Clindamycin	≥1	R
Benzylpenicillin (Oral)	≥8	R	Erythromycin	≥8	R
Benzylpenicillin (Other)	≥8	N/A	Levofloxacin	0.5	S
Benzylpenicillin (Pneumonia)	≥8	R	Linezolid	≤2	S
Cefotaxime (Meningitis)	4	R	Moxifloxacin	0.12	S
Cefotaxime (Other)	4	R	Tetracycline	≥16	R
Ceftriaxone (Meningitis)	4	R	Tigecycline	≤0.06	S
Ceftriaxone (Other)	4	R	Trimethoprim/Sulfamethoxazole	160	R
			Vancomycin	0.5	S

^a Antibiotic susceptibility was obtained using a Vitek2 AST-ST02 cards.


^b Parameter set: MIC Interpretation Guideline: CLSI M100-S27 (2017)

Therapeutic Interpretation Guideline: Natural Resistance

ATCC® Genome Portal – Ways to view AMR

View / Download Phenotypical Data



HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION LOG IN Become a Supporting Member

Streptococcus > **Streptococcus pneumoniae**

Streptococcus pneumoniae (ATCC® BAA-3330™)

Overview Genome Browser Related Genomes Quality Control

Sequencing Quality Control

Sequencing Quality Control Metrics	
Number of trimmed reads	1,908,100
Median Q score, all bases	40
Percent of median Q scores per cycle greater than 25	1%
Ambiguous content (% N bases)	0%

Assembly Quality Control

Assembly Quality Control Metrics	
Estimated genome completeness	99.62%
Estimated genome contamination	0.31%
Illumina depth of coverage	238.621x
ONT depth of coverage	864.051x



Notes

This bacterium was assembled, annotated, and curated using Oatmeal v3.0. See our technical documentation for further details: <https://docs.onecodex.com/en/articles/5357065-assembly-qc-and-annotation-pipeline-versions>. This item includes additional phenotypic antimicrobial resistance (AMR) data. You may download the Vitek2 antibiotic susceptibility report here: <https://www.atcc.org/-/media/product-assets/documents/antimicrobial-susceptibility-reports/BAA-3330-antibiotic-susceptibility-profile.pdf>

ATCC® Genome Portal – Ways to view AMR

View / Download Phenotypical Data





[Applications](#) [Cell Products](#) [Microbe Products](#) [Services](#) [Federal Solutions](#) [Resources](#) [About Us](#)

[Home](#) > [Microbe Products](#) > [Bacteriology and Archaea](#) > [BAA-3330](#)

Streptococcus pneumoniae (Klein) Chester

BAA-3330™

[Download Genome](#)
[LEARN ABOUT THE ATCC GENOME PORTAL >](#)

This drug-resistant bacterium was isolated in 2013 from the respiratory bronchoalveolar lavage of a 22-year-old male in Italy. This product can be used for all stages of the discovery and development process for novel antimicrobials and therapeutics, molecular-based detection assays, and updated sterility protocols.

Product category	Bacteria
Product type	Drug-resistant bacterium
Strain designation	978059
Type strain	No
Genome sequenced strain	Yes

Buy Now

Price: \$304.00 ea

Discounts may be available for our fellow nonprofit organizations. [Login](#) to see your price.

Generally ships within 1-3 business days

Quantity

[Add to cart](#)

[Add to list](#)

Documentation

[Product sheet](#) [Certificate of analysis](#) [Safety data sheet](#)

[Antimicrobial resistance report](#) [Characterization data](#)

ATCC® Genome Portal – Ways to view AMR

Download Resistome Dataset - GUI



The screenshot shows the ATCC Genome Portal interface for *Streptococcus pneumoniae* (ATCC® BAA-3330™). The page includes a navigation bar with links for HOME, MICROBE PRODUCTS, CELL PRODUCTS, and DOCUMENTATION, along with a search bar and a user profile dropdown for JPETRONE@ATCC.ORG. The main content area features a breadcrumb trail, the organism name, and tabs for Overview, Genome Browser, Related Genomes, and Quality Control. A 'DOWNLOAD DATA' button is highlighted, opening a dropdown menu with four data options:

- Assembly Data** (1 dataset available) - 592.8 kB
File: *Streptococcus_pneumoniae_ATCC_BAA_3330.fasta*
Date: October 31, 2025
Description: A FASTA file containing assembled genomic sequences
- Annotations Data** (1 dataset available) - 1.5 MB
File: *Streptococcus_pneumoniae_ATCC_BAA_3330.gbk*
Date: October 31, 2025
Description: A Genbank (.gbk) file containing annotations for the associated assembly
- Methylation Data** (1 dataset available) - 4.5 MB
File: *BAA-3330_f30b5b23cc3d401a88ebc98588f2ab1d_methyl_beds.zip*
Date: July 11, 2025
Description: A JSON file containing metadata for the Methylation experiments and the associated data in BED format
- Resistome Data** (1 dataset available) - 619.8 kB
File: *BAA-3330_resistome_dataset.zip*
Date: November 4, 2025
Description: A spreadsheet of antimicrobial resistance phenotypes

A 'Download' button is located at the bottom right of the dropdown menu. To the right of the dropdown, an 'Assembly Summary' table provides key metrics:

Assembly Summary	
Date Published	June 11, 2025
Length	2,130,590 nt
Sequencing Technology	Illumina + Oxford Nanopore Hybrid Assembly
Number of Contigs	1 (All Circularized)
N50	2,130,590 nt
%GC	39.85%

Below the assembly summary is an 'Annotations Summary' table:

Annotations Summary	
Number of CDS	2148
Number of Hypothetical Proteins	217
Number of tRNA	59
Number of 5s rRNA	4
Number of 16s rRNA	4
Number of 23s rRNA	4

ATCC[®] Genome Portal – Ways to view AMR

Download Resistome Dataset - GUI



1	Locus Tag	RESFINDER HIT	AMRFINDER HIT	CARD HIT	PGAP Gene	PGAP Sequence Name	AMRFinder Class	CARD Sequence Name	CARD Resistance Mechanism	
1783	ATCC_BAA3330_001821	True	True	True	erm(B)	23S rRNA (adenine(2058)-N(6))-methyltransferase Erm(B)	LINCOSAMIDE/MACROLIDE/STREPTOGRAMIN	ErmB	antibiotic target alteration	erythromycin; roxithromycin; lincomycin; telithromycin; clarithromycin
1784	ATCC_BAA3330_001822					ISL3 family transposase				
1785	ATCC_BAA3330_001823					hypothetical protein				
1786	ATCC_BAA3330_001824					antirestriction protein ArdA				
1787	ATCC_BAA3330_001825					conjugal transfer protein				
1788	ATCC_BAA3330_001826					ATP-binding protein				
1789	ATCC_BAA3330_001827					CD3337/EF1877 family mobilome membrane protein				
1790	ATCC_BAA3330_001828					lysozyme family protein				
1791	ATCC_BAA3330_001829					conjugal transfer protein				
1792	ATCC_BAA3330_001830					hypothetical protein				
1793	ATCC_BAA3330_001831	True	True	True	tet(M)	tetracycline resistance ribosomal protection protein Tet(M)	TETRACYCLINE	tet(M)	antibiotic target protection	tetracycline; doxycycline; minocycline; chlortetracycline; demeclocycline
1794	ATCC_BAA3330_001832					hypothetical protein				
1795	ATCC_BAA3330_001833					imsC-like protein				
1796	ATCC_BAA3330_001834					hypothetical protein				
1797	ATCC_BAA3330_001835					DUF5960 family protein				
1798	ATCC_BAA3330_001836	True	True		msr(D)	ABC-F type ribosomal protection protein Msr(D)	MACROLIDE/STREPTOGRAMIN			
1799	ATCC_BAA3330_001837	True	True		mef(A)	macrolide efflux MFS transporter Mef(A)	MACROLIDE			
1800	ATCC_BAA3330_001838					DNA recombinase				
1801	ATCC_BAA3330_001839					cysteine-rich KTR domain-containing protein				
1802	ATCC_BAA3330_001840					helix-turn-helix domain-containing protein				
1803	ATCC_BAA3330_001841					RNA polymerase sigma factor				
1804	ATCC_BAA3330_001842					helix-turn-helix domain-containing protein				
1805	ATCC_BAA3330_001843					excisionase				
1806	ATCC_BAA3330_001844					tyrosine-type recombinase/integrase				
1807	ATCC_BAA3330_001845					DUF1430 domain-containing protein				
1808	ATCC_BAA3330_001846					ABC transporter ATP-binding protein				
1809	ATCC_BAA3330_001847				ndk	nucleoside-diphosphate kinase				
1810	ATCC_BAA3330_001848				rpoC	DNA-directed RNA polymerase subunit beta'				
1974	ATCC_BAA3330_002028					LysM peptidoglycan-binding domain-containing protein				
1975	ATCC_BAA3330_002029					HAD family hydrolase				
1976	ATCC_BAA3330_002030					MATE family efflux transporter				
1977	ATCC_BAA3330_002031				thrC	threonine synthase				
1978	ATCC_BAA3330_002032					gamma-glutamylcysteine synthetase				
1979	ATCC_BAA3330_002049					cytidine deaminase family protein				
1980	ATCC_BAA3330_002050				gltX	glutamate--trNA ligase				

```
{
  "amr_software_metadata": {
    "annotation_software_version": "version 2025-05-06.build7983",
    "AMRFinderPlus_version": "4.0.23",
    "AMRFinderPlus_database_version": "2025-07-16.1",
    "ResFinder_version": "4.7.2",
    "RGI_version": "6.0.4",
    "CARD_database_version": "4.0.1",
    "date_annotations_pan": "09/03/2025 01:08:39",
    "annotation_software": "PGAP"
  },
  "product_metadata": {
    "ATCC_catalog_number": "BAA-3330",
    "ATCC_item_name": "Streptococcus pneumoniae",
    "lot_number": "700"
  },
  "assembly_metadata": {
    "assembly_id": "939abb8ff2534c01",
    "genome_id": "391623df94964877"
  }
}
```

Linking Phenotype and Genotype of ATCC Strains



Antibiotic Susceptibility of *Streptococcus pneumoniae* Strain 978059 ATCC® BAA-3330™

Antimicrobial	MIC ^a	Interpretation ^b	Antimicrobial	MIC	Interpretation
Benzylpenicillin (Meningitis)	≥8	R	Clindamycin	≥1	R
Benzylpenicillin (Oral)	≥8	R	Erythromycin	≥8	R
Benzylpenicillin (Other)	≥8	N/A	Levofloxacin	0.5	S
Benzylpenicillin (Pneumonia)	≥8	R	Linezolid	≤2	S
Cefotaxime (Meningitis)	4	R	Moxifloxacin	0.12	S
Cefotaxime (Other)	4	R	Tetracycline	≥16	R
Ceftriaxone (Meningitis)	4	R	Tigecycline	≤0.06	S
Ceftriaxone (Other)	4	R	Trimethoprim/Sulfamethoxazole	160	R
			Vancomycin	0.5	S

← VITEK Phenotype

^a Antibiotic susceptibility was obtained using a Vitek2 AST-ST02 cards.

^b Parameter set: MIC Interpretation Guideline: CLSI M100-S27 (2017)
Therapeutic Interpretation Guideline: Natural Resistance

Resistome Genotype ↓

Locus Tag	RESFINDER HIT	AMRFINDER HIT	CARD HIT	PGAP Gene symbol	AMRFinder Subclass	CARD Sequence Name	CARD Resistance Mechanism	CARD Antibiotic	RESFINDER Phenotype
ATCC_BAA3330_000368		True		pbp2X	BETA-LACTAM				
ATCC_BAA3330_000624			True	ldcB		vanY gene in vanB cluster	antibiotic target alteration	vancomycin	
ATCC_BAA3330_001083		True		aac(6')	AMINOGLYCOSIDE				
ATCC_BAA3330_001167	True	True	True		QUINOLONE	pmrA	antibiotic efflux	ciprofloxacin; norfloxacin	
ATCC_BAA3330_001533		True		pbp2b	BETA-LACTAM				
ATCC_BAA3330_001821	True	True	True	erm(B)	CLINDAMYCIN/ERYTHROMYCIN/STREP TOGRAMIN B	ErmB	antibiotic target alteration	erythromycin; roxithromycin; lincomycin; telithromycin; clarithromycin; clindamycin; tylosin; spiramycin; azithromycin; dirithromycin; pristinamycin IA; quinupristin; virginiamycin M1; madumycin II; griseoviridin; dalfopristin; pristinamycin IB; virginiamycin S2; pristinamycin IC; vernamycin C; patricin A; patricin B; ostreogrycin B3; oleandomycin	Erythromycin, Lincomycin, Clindamycin, Quinupristin, Pristinamycin IA, Virginiamycin S
ATCC_BAA3330_001831	True	True	True	tet(M)	TETRACYCLINE	tet(M)	antibiotic target protection	tetracycline; doxycycline; minocycline; chlortetracycline; demeclocycline; oxytetracycline	Doxycycline, Tetracycline, Minocycline
ATCC_BAA3330_001836	True	True		msr(D)	AZITHROMYCIN/ERYTHROMYCIN/STREPTOGRAMIN B				Erythromycin, Azithromycin, Telithromycin, Quinupristin, Pristinamycin IA, Virginiamycin S
ATCC_BAA3330_001837	True	True		mef(A)	AZITHROMYCIN/ERYTHROMYCIN				Erythromycin, Azithromycin
ATCC_BAA3330_002054			True	patB		patB	antibiotic efflux	ciprofloxacin; norfloxacin	
ATCC_BAA3330_002056			True	patA		patA	antibiotic efflux	ciprofloxacin; norfloxacin	
ATCC_BAA3330_002083			True			RlmA(II)	antibiotic target alteration	tylosin; mycinamicin	

ATCC® Genome Portal – Download with Lot Number (All Users)



ATCC®

HOME MICROBE PRODUCTS CELL PRODUCTS DOCUMENTATION Quick Search...

Welcome to the ATCC Genome Portal

The only authenticated reference omics database for ATCC biom...

VIEW ALL MICROBES >

VIEW ALL CELL LINES >

Materials and Methods
REST API
How to Cite the ATCC Genome Portal
Download Genome for Purchased Strain
Release Notes
Explore Our Resources

*** > Nakaseomyces > *Candida glabrata*

Candida glabrata (ATCC® 2001™)

Overview Annotations Related Genomes Quality Control

DOWNLOAD ASSEMBLY DOWNLOAD ANNOTATIONS

Assembly Summary

Date Published	October 30, 2020
Length	12,583,835 nt
Sequencing Technology	Illumina + Oxford Nanopore
Number of Contigs	20
Assembly Level	Scaffold
N50	1,108,508 nt
%GC	38.77%

Annotations Summary Input Reads Summary

Details required to download data


Downloads are only available to ATCC Genome Portal supporting members or those who have purchased a corresponding physical product.

You can select a supporting membership plan or enter the lot number associated with your ATCC product below to continue.

Lot Number

Please enter lot number

Cancel View Plans



Access to the entire database...purchase an annual *Supporting Membership*



	Free	Individual	Research Group	Institution
View organism and genome metadata, assemblies, and annotations	✓	✓	✓	✓
Search for genomes of interest	✓	✓	✓	✓
Purchase the corresponding authenticated ATCC source materials	✓	✓	✓	✓
Download genome assemblies and annotations	Only for purchased products	All products	All products	All products
Access our REST-API	Not available	✓	✓	✓
Analyze isolates with Discrepancy Reports	Fee for each report	12 free reports per year	60 free reports per year	Inquire
Members with full access	0	1	5	Unlimited
		\$600/\$1,800 NFP / Profit	\$2,400/\$7,200 NFP / Profit	Inquire

Meet ATCC Scientists at ASM Microbe

Discover trusted microbial standards, reference genomes, and expert guidance for reproducible results



Visit us at
booth #1724



Check out our 10
scientific posters



Explore our
resources



See our full agenda for ASM Microbe 2026