

ATCC Medium: 1256 Sulfolobus Medium

(NH ₄) ₂ SO ₄	1.3 g
KH ₂ PO ₄	0.28 g
MgSO ₄ x 7H ₂ O.....	0.25 g
CaCl ₂ x 2H ₂ O	0.07 g
FeCl ₃ x 6H ₂ O.....	0.02 g
MnCl ₂ x 4H ₂ O.....	1.8 mg
Na ₂ B ₄ O ₇ x 10H ₂ O.....	4.5 mg
ZnSO ₄ x 7H ₂ O	0.22 mg
CuCl ₂ x 2H ₂ O.....	0.05 mg
Na ₂ MoO ₄ x 2H ₂ O.....	0.03 mg
VOSO ₄ x 2H ₂ O.....	0.03 mg
CoSO ₄	0.01 mg
Yeast Extract.....	1.0 g
DI Water.....	1000 ml

Adjust pH to 2.0 with 10 N H₂SO₄ (at room temperature). Filter-sterilize and dispense into appropriate vessel.