

**ATCC medium: 2455 Middlebrook 7H9 with OADC-Pyruvate-Tryptone**

Middlebrook 7H9 Broth (BD 271310).....4.7 g  
Glycerol.....2.0 ml  
Tryptone.....1.0 g  
Sodium pyruvate.....4.0 g  
OADC Enrichment (BD 211886).....100.0 ml  
Agar (if needed).....15.0 g  
Distilled deionized water.....900.0 ml

Combine all ingredients except the OADC Enrichment. Mix well. Autoclave for 10 minutes at 121C (see note). Aseptically add OADC Enrichment to the cooled basal medium and dispense as required. If agar is desired, add 15 grams per liter to the basal medium and boil before autoclaving.

Note: Difco Manual, 11<sup>th</sup> ed., calls for 10 minutes at 121C. Be sure to autoclave for the stated time.