

ATCC Medium: 2832 Reduced YPD Medium

KH ₂ PO ₄	1.0 g
MgSO ₄ x 7H ₂ O.....	0.5 g
Casein peptone.....	5.0 g
Yeast Extract.....	1.0 g
Glucose.....	10.0 g
Wheat germ oil (Sigma).....	3.0 g
Cysteine x HCl.....	0.3 g
DI Water.....	1000 ml

Adjust pH to 6.8-7.2. Fill tubes to the top (approx.. 15.5 ml/tube). Autoclave at 121C.