

ATCC medium: 178 *Pseudomonas* medium

K₂HPO₄1.0 g
CaCl₂0.5 g
MgSO₄ . 7H₂O0.5 g
Na₂SO₄0.25 g
MnSO₄ . 4H₂O0.05 g
FeSO₄ . 7H₂O0.05 g
Na₂MoO₄ . 2H₂O.....0.05 g
Mannitol.....10.0 g
Agar.....16.0 g
Distilled water.....1.0 L

Adjust pH to 7.4-7.6
Autoclave at 121C for 15 minutes.