

**ATCC: Medium 2694 PYE medium for 43427**

Yeast Extract.....	10.0g
Sodium Pyruvate.....	1.1g
K <sub>2</sub> HPO <sub>4</sub> .....	1.0g
MgSO <sub>4</sub> ·7H <sub>2</sub> O.....	0.5g
Agar (if required).....	20.0g
Distilled Water .....	1.0L

Adjust pH to 7.0 – 7.2

Dispense into desired vessel and autoclave @ 121°C.

\*\*For solid medium add agar to a final concentration of 2%.