

**ATCC medium: 1244 Van Niel's yeast agar (ATCC medium 112) with 5 mM glutamate**

ATCC Medium 112 (see below) with 5 mM glutamate

*ATCC Medium 112:*

K <sub>2</sub> HPO <sub>4</sub> .....	1.0 g
MgSO <sub>4</sub> .....	0.5 g
Yeast extract.....	10.0 g
Agar.....	20.0 g
Tap water.....	1.0 L

Adjust pH to 7.0-7.2. Autoclave at 121C for 15 minutes.