

ATCC medium: 1160 Dubos agar with filter paper

NaNO₃0.5 g
K₂HPO₄1.0 g
MgSO₄ . 7H₂O0.5 g
KCl.....0.5 g
FeSO₄ . 7H₂O0.01 g
Agar, Bacto (BD 214050).....15.0 g
Distilled water.....1.0 L

Adjust pH to 7.2 after autoclaving. Sterile strips of Whatman No.1 filter paper are laid on surface of agar. A strip of filter paper is inserted in each tube of broth.