

ATCC medium: 2003 Czapek yeast extract agar

K₂HPO₄1.0 g
Czapek Concentrate (see below)....10.0 ml
Yeast extract.....5.0 g
Sucrose.....30.0 g
Agar.....15.0 g
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.

Czapek Concentrate:

NaNO₃30.0 g
KCl.....5.0 g
MgSO₄ . 7H₂O5.0 g
FeSO₄ . 7H₂O0.1 g
ZnSO₄ . 7H₂O0.1 g
CuSO₄ . 5H₂O0.05 g
Distilled water.....1.0 L