

ATCC medium: 711 PYB

Proteose Peptone (BD 211684).....1.0 g
Yeast extract.....1.0 g
Agar (if needed).....20.0 g
Distilled water.....900.0 ml

Prepare and sterilize separately each of the following stock solutions and add to the basal medium as indicated below to avoid precipitation. Make final volume to 1.0 L.

Stock Solutions:

0.05 M CaCl₂4.0 ml
0.4 M MgSO₄ · 7H₂O2.5 ml
0.25 M Na₂HPO₄8.0 ml
0.25 M KH₂PO₄32.0 ml

pH should be 6.5. Autoclave 25 minutes at 121C.