

**ATCC medium: 668 SMC medium**

PPL0 Broth w/o CV (Mycoplasma Broth) (BD 255420).....21.0 g  
Tryptone (BD 211705).....10.0 g  
Sorbitol.....70.0 g  
Distilled water.....650.0 ml

Adjust pH to 7.5.  
Autoclave at 121C for 15 minutes.

Filter-sterilize and add aseptically:

Horse serum.....200.0 ml  
Yeast Extract Solution (GIBCO 360-8180).....100.0 ml  
50% Fructose solution.....2.0 ml  
50% Sucrose solution.....20.0 ml  
50% Glucose solution.....2.0 ml  
42% L-Arginine HCl.....10.0 ml  
0.1% Phenol red solution.....20.0 ml