

**ATCC medium: 1381 PPYG medium**

Peptone.....5.0 g  
Yeast extract.....1.5 g  
Glucose.....5.0 g  
NaCl.....1.5 g  
Na<sub>2</sub>HPO<sub>4</sub> · 12H<sub>2</sub>O.....1.5 g  
MgCl<sub>2</sub> · 6H<sub>2</sub>O .....0.1 g  
Na<sub>2</sub>CO<sub>3</sub> .....5.03 g  
Agar (if needed).....15.0 g  
Distilled water.....1.0 L

Sterilize glucose and Na<sub>2</sub>CO<sub>3</sub> separately and adjust final pH to 10.5-11.0. Autoclave at 121C for 15 minutes.