

ATCC medium: 2169 *Halobaculum gomorrense* medium

NaCl.....125.0 g
MgCl₂ . 6H₂O160.0 g
K₂SO₄5.0 g
CaCl₂ . 2H₂O0.1 g
Yeast extract.....1.0 g
Casamino acids.....1.0 g
Starch.....2.0 g
Distilled water.....1.0 L

Adjust pH to 7.0.

Autoclave at 121C for 15 minutes.

Note:

If agar is called for, the high salt concentration may require that the agar and other ingredients be autoclaved separately.