

ATCC medium: 1421 LHET2 medium (ATCC medium 1168) with 0.01% yeast extract or yeast autolysate

ATCC Medium 1168 (see below) with 0.01% yeast extract or yeast autolysate

ATCC Medium 1168:

Solution A:

(NH ₄) ₂ SO ₄	2.0 g
KCl.....	0.1 g
K ₂ HPO ₄	0.5 g
MgSO ₄ . 7H ₂ O	0.5 g
Trypticase Soy Broth without Dextrose (BD 211774).....	0.1 g
Distilled water.....	500.0 ml

Adjust pH to 2.5-3.0 using 1 N H₂SO₄.

Solution B:

Glucose.....	1.0 g
Agar (for solid medium).....	12.0 g
Distilled water.....	500.0 ml

Autoclave Solutions A and B separately. Cool to 50C. Add Solution A to Solution B.