

ATCC medium: 2676 Natrinema Medium

Yeast Extract.....5.00 g
Casamino acids.....5.00 g
Na-glutamate.....1.00 g
NaCl.....200.00 g
KCl.....2.00 g
MnCl₂ x 4 H₂O.....5.00 g
FeCl₂ x 4 H₂O.....5.00 g
MgSO₄ x 7 H₂O.....20.00 g
Na₃ citrate.....3.00 g
Agar.....15.00 g
Distilled water.....1.00 L

Adjust pH to 7.0 - 7.2 Autoclave at 121C for 15 minutes.