

ATCC medium: 566 *Megasphaera* medium

KH ₂ PO ₄	1.6 g
K ₂ HPO ₄	3.2 g
Sodium lactate (60% solution).....	16.0 ml
Yeast extract.....	4.0 g
NH ₄ Cl	0.5 g
CaCl ₂	0.2 g
MgCl ₂	0.2 g
Sodium thioglycollate.....	0.45 g
Agar.....	1.0 g
Distilled water to.....	1.0 L

Adjust pH to 7.0.

Autoclave at 121C for 15 minutes.