

ATCC Medium: 2847 GYP w/ salts

Base Medium

Yeast Extract.....	10 g
Peptone.....	10 g
Sodium acetate.....	10 g
Salt B Solution (see below).....	5 ml
DI Water.....	895 ml

Adjust pH to 6.8. Autoclave at 121C. When medium is cooled, aseptically add the following filter sterilized solution:

Glucose Solution

Glucose.....	20 g
DI Water.....	100 ml

Salt B Solution

MgSO ₄ x 7 H ₂ O.....	4 g
FeSO ₄ x 7H ₂ O.....	0.2 g
NaCl.....	0.2 g
MnSO ₄ x 4H ₂ O.....	0.16 g
DI Water.....	100 ml