

**ATCC medium: 666 Modified Wakimoto medium**

Ca(NO<sub>3</sub>)<sub>2</sub> · 4H<sub>2</sub>O.....0.5 g  
Na<sub>2</sub>HPO<sub>4</sub> · 12H<sub>2</sub>O.....2.0 g  
Peptone.....5.0 g  
Sucrose.....15.0 g  
FeSO<sub>4</sub> · 7H<sub>2</sub>O .....0.5 g  
Agar.....15.0 g  
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.