

ATCC Medium 2441: King's Medium B w/ 50 ug/mL Rifampicin

Proteose Peptone No. 3.....	20.0 g
Glycerol.....	10.0 mL
K ₂ HPO ₄	1.5 g
MgSO ₄ x 7H ₂ O.....	1.5 g
Agar (if required).....	15.0 g
DI Water.....	950.0 mL
Rifampicin Solution (See below).....	50.0 mL

Mix well and adjust to pH 7.2 +/- 0.2. Autoclave at 121° C for 15 minutes. Cool to 55° C then add sterile Rifampicin solution.

Rifampicin Solution

Rifampicin.....	50.0 mg
DI Water.....	50.0 mL

Filter sterilize.