

**ATCC medium: 1287 Kerosene-mineral salts medium (Bushnell-Haas)**

MgSO<sub>4</sub> . 7H<sub>2</sub>O .....0.2 g  
CaCl<sub>2</sub> .....0.02 g  
KH<sub>2</sub>PO<sub>4</sub> .....1.0 g  
NH<sub>4</sub>NO<sub>3</sub> .....1.0 g  
K<sub>2</sub>HPO<sub>4</sub> .....1.0 g  
FeCl<sub>3</sub> concentrated solution (15 g/25 ml distilled water) ...2 drops  
Distilled water to.....1.0 L

Adjust pH to 6.9 - 7.0 with dilute NaOH before sterilization.  
Overlay sterile kerosene on medium at time of use (2%).