

ATCC medium: 2314 R2A medium with NaCl

Prepare ATCC Medium 2258 (see below); add 20.0 g of NaCl to the other components listed.

ATCC Medium 2258:

Yeast extract.....	0.5 g
Proteose Peptone No. 3 (BD 211693).....	0.5 g
Casamino acids.....	0.5 g
Glucose.....	0.5 g
Soluble starch.....	0.5 g
Sodium pyruvate.....	0.3 g
K ₂ HPO ₄	0.3 g
MgSO ₄ . 7H ₂ O	0.05 g
Agar.....	15.0 g
Distilled water.....	1.0 L

Adjust final pH to 7.2 with potassium phosphate before adding the agar.
Autoclave at 121C for 15 minutes.