

ATCC Medium 671: Dictyostelium Medium

Peptone.....	14.3 g
Yeast Extract (oxid).....	7.15 g
Na ₂ HPO ₄ x 12H ₂ O.....	1.28 g
KH ₂ PO ₄	0.49 g
Agar.....	15.0 g
DI Water.....	950.0 mL
Glucose Solution (See Below).....	50.0 mL

Adjust pH to 6.7

Mix and Autoclave at 121°C.

Glucose Solution

Glucose (Dextrose).....	15.4 g
DI Water.....	50.0 mL

Filter sterilize. Aseptically add to base media.