

Antibiotic Susceptibility of *Pseudomonas aeruginosa* ATCC® BAA-2111™

Antimicrobial	MIC ^a	Interpretation ^b	Antimicrobial	MIC	Interpretation
Amikacin	8	S	Imipenem	0.5	S
Aztreonam	2	S	Levofloxacin	≤1	S
Cefepime	≤4	S	Meropenem	≤0.5	S
Ceftazidime	1	S	Piperacillin	≤16	S
Ciprofloxacin	≤0.5	S	Ticarcillin / clavulanic acid	16/2	S
Doripenem	≤0.5	S	Tobramycin	≤2	S
Gentamicin	2	S			

^a Antibiotic susceptibility was obtained using the Sensititre® brand of antimicrobial titration analysis system.

^b Parameter set: MIC Interpretation Guideline: CLSI M100 ED32:2022