

**ATCC medium: 1106 PPLO broth with bovine serum**

PPLO Broth Base w/o CV (Mycoplasma Broth Base) (BD 211458).....21.0 g  
Phenol red solution, 1% stock solution.....2.0 ml  
Deionized water.....865.0 ml

Adjust pH to 7.5. Autoclave.

Aseptically add:

Glucose (20% solution)\*.....25.0 ml  
Bovine serum\*.....10.0 ml  
Yeast Extract Solution (GIBCO 18180).....100.0 ml

\*Filter-sterilize.