

**ATCC medium: 1108 *Renibacterium* KDM-2 medium**

Peptone.....1.0 % w/v  
Yeast extract.....0.05 % w/v  
L-Cysteine . HCl.....0.1 % w/v  
Agar.....1.5 % w/v  
Fetal bovine serum....20.0 % w/v

Dissolve peptone, yeast extract and cysteine in 80% final volume distilled water. Adjust to pH 6.5 with NaOH. Add agar and dissolve by heating. Autoclave at 121C for 15 minutes. Cool to 45C, add serum at 45C.