

ATCC medium: 1143 Yeast synthetic medium with 5'-dTMP

Yeast Nitrogen Base w/o amino acids (BD 291940).....6.7 g
Casamino Acids, Vitamin Assay (BD 228830).....2.0 g
D (+) glucose.....20.0 g
2'-Deoxythymidine-5'-monophosphate.....7.2 mg
Agar.....15.0 g
Distilled water.....1.0 L

Dissolve first four components in 200 ml distilled water and filter-sterilize. Dissolve and autoclave agar in 800 ml distilled water. Combine both solutions aseptically.