

ATCC medium: 1144 Xylose-YP medium

Xylose.....	10.0 g
Yeast Extract (BD 212750).....	10.0 g
Peptone, Bacto (BD 211677).....	10.0 g
MgSO ₄ . 7H ₂ O	200.0 mg
MnSO ₄ . 4H ₂ O	10.0 mg
FeSO ₄ . 7H ₂ O	10.0 mg
NaCl.....	10.0 mg
Distilled water.....	1.0 L

Adjust final pH to 6.8.

Note: Add 20 g/L of CaCO₃ when preparing agar plates. Resuspend CaCO₃ and pour agar immediately after autoclaving.