

**ATCC medium: 1216 Modified SMC**

PPL0 broth w/o CV (BD 255420).....21.0 g  
Tryptone.....10.0 g  
Sorbitol.....70.0 g  
NaOH, 1.0 M.....6.0 ml  
Distilled water.....870.0 ml

Autoclave at 121C for 15 minutes. Cool and aseptically add the following filter-sterilized solutions:

*Solution 1:*

Glucose.....1.0 g  
Fructose.....1.0 g  
Sucrose.....10.0 g  
TC Yeastolate (BD 255772).....2.0 g  
Phenol red.....20.0 mg  
Distilled water.....100.0 ml

*Solution 2:*

1% Bovine serum albumin fraction V.....10.0 ml

*Solution 3:*

Horse serum (heat-inactivated).....20.0 ml

Adjust medium for final pH 7.7-7.8.