

**ATCC Medium: 1271 Benzoate Medium**

(NH <sub>4</sub> ) <sub>2</sub> HPO <sub>4</sub> .....	3.0 g
KH <sub>2</sub> PO <sub>4</sub> .....	1.2 g
NaCl.....	5.0 g
MgSO <sub>4</sub> x 7H <sub>2</sub> O.....	0.2 g
Yeast Extract.....	0.5 g
Agar, Noble (BD 214230).....	20.0 g
DI Water.....	950 ml

Combine all ingredients except sodium benzoate. Autoclave at 121°C. Aseptically add the following filter-sterilized sodium benzoate solution:

**Sodium Benzoate Solution**

Sodium benzoate.....	3.0 g
DI Water.....	50 ml