

ATCC medium: 1513 PY-CMC medium

Carboxymethyl cellulose.....10.0 g
Polypeptone Peptone (BD 211910).....5.0 g
Yeast extract.....5.0 g
NaCl.....5.0 g
KH₂PO₄1.0 g
MgSO₄ . 7H₂O2.0 g
Na₂CO₃*10.0 g
Agar (if necessary).....15.0 g
Distilled water.....1.0 L

Adjust medium, if necessary, for final pH 9.5. Autoclave at 121C for 15 minutes.

*Sterilize sodium carbonate as a separate solution and add aseptically to the remainder of the medium.