

ATCC Medium: 1611 U4 Medium

Hanks' Balanced Salt Solution, 10X (see below).....	4.0 ml
Hartley Digest Broth (HIMEDIA M551).....	20.0 ml
Fetal Bovine Serum.....	15.0 ml
Yeast Extract Solution (GIBCO 18180).....	10.0 ml
1% Phenol red.....	0.2 ml
20% Urea.....	0.25 ml
0.025% MgSO ₄ x 5H ₂ O.....	1.0 ml
DI Water to.....	100.0 ml

Adjust pH to 6.0-6.2 with HCl and filter-sterilize complete medium.

Hanks' Balanced Salt Solution, 10X

CaCl ₂	1.4 g
KCl.....	4.0 g
KH ₂ PO ₄	0.6 g
MgCl ₂ x 6H ₂ O	1.0 g
MgSO ₄ x 7H ₂ O	1.0 g
NaCl.....	80.0 g
Na ₂ HPO ₄ x 7H ₂ O.....	0.9 g
Glucose.....	10.0 g
DI Water to.....	1000.0 ml