

**ATCC medium: 1674 GC agar base with Kellogg's supplement**

GC Agar Base (BD 212307).....36.0 g  
Kellogg's Supplement (see below).....10.0 ml  
Distilled water.....990.0 ml

Autoclave basal medium at 121C for 15 minutes; cool and add supplement aseptically.

*Kellogg's Supplement:*

Glucose.....40.0 g  
Glutamine.....0.5 g  
Ferric nitrate.....0.05 g  
Coccarboxylase.....2.0 mg  
Distilled water.....100.0 ml

Filter-sterilize.