

ATCC medium: 1676 Modified Van Niel's medium

MgSO₄0.1 g
EDTA.....2.0 mg
Yeast extract.....10.0 g
Trace Elements (see below).....10.0 ml
4% K₂HPO₄ (add aseptically after sterilization)2.5 ml
Distilled water.....985.0 ml

Adjust medium for final pH 7.1. Autoclave at 121C for 15 minutes.

Trace Elements:

Ferric ammonium citrate.....0.2 g
CaCl₂ . 2H₂O0.3 g
Distilled water.....100.0 ml

Filter-sterilize and add aseptically to basal medium.