

ATCC medium: 1678 YGCP medium

Yeast extract.....	2.5 g
Glucose.....	20.0 g
Peptone, Bacto (BD 211677).....	2.5 g
NaCl.....	1.0 g
K ₂ HPO ₄	1.0 g
MgSO ₄ (anhydrous)	0.5 g
CaCO ₃	10.0 g
Agar.....	17.5 g
Distilled water.....	1.0 L

Add components except CaCO₃ and adjust to pH 7.2. Add CaCO₃ and autoclave at 121C for 15 minutes.