

**ATCC medium: 1687 *Flavobacterium* medium**

K<sub>2</sub>HPO<sub>4</sub> .....0.65 g  
KH<sub>2</sub>PO<sub>4</sub> .....0.19 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....0.1 g  
NaNO<sub>3</sub> .....0.5 g  
Sodium glutamate.....4.0 g  
0.01 M FeSO<sub>4</sub>\* .....2.0 ml  
Distilled water.....1.0 L

Adjust medium for final pH 7.4. Autoclave at 121C for 15 minutes.

\*Add aseptically after autoclaving.