

**ATCC medium: 1763 YM-catalase agar**

Yeast extract.....5.0 g  
Malt extract.....5.0 g  
K<sub>2</sub>HPO<sub>4</sub> .....5.74 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....205.0 mg  
NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub>.....1.15 g  
Catalase.....60.0 mg  
Agar.....15.0 g  
Distilled water.....1.0 L

Filter-sterilize catalase and MgSO<sub>4</sub> separately. Autoclave remainder of medium at 121C for 15 minutes.