

1877

ISP #1 Medium

ISP#1 (BD 276910)	8 g
DI Water	1000 mL

Autoclave at 121 C for 15 minutes.

To make agar:

Tryptone	5.0	g
Yeast extract	3.0	g
Agar	15.0	g
Distilled water	1000	mL

Adjust medium for final pH 7.0 – 7.2.
Autoclave at 121 C for 15 minutes.

To make from scratch:

Tryptone	5.0	g
Yeast extract	3.0	g
Distilled water	1000	mL

Adjust medium for final pH 7.0 – 7.2.
Autoclave at 121 C for 15 minutes.