

ATCC medium: 1906 Celery decoction agar

Celery Decoction (see below).....320.0 ml
Potato Dextrose Broth (BD 254920).....4.0 g
Agar.....15.0 g
Distilled water.....680.0 ml

Adjust to pH 5.6 with 20% H₂SO₄. Autoclave at 121C for 15 minutes.

Celery Decoction:

Celery leaves and petioles.....50.0 g
Distilled water.....500.0 ml

Coarsely cut celery into pieces, combine with water and puree in blender at high speed for 2 minutes. Strain through a kitchen-type sieve. Boil the mixture for 5 minutes and strain through a very fine sieve or muslin.