

ATCC medium: 1949 LB medium with chloramphenicol and ampicillin

ATCC Medium 1065 (see below) with 50.0 mcg/ml chloramphenicol and 100.0 mcg/ml ampicillin

Filter-sterilize solutions of chloramphenicol and ampicillin and add to sterile, cooled medium to achieve the stated final concentration.

ATCC Medium 1065:

Tryptone (BD 211705).....	10.0 g
Yeast Extract (BD 212750).....	5.0 g
NaCl.....	10.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.