

ATCC Medium 196: Yeast Malt Extract Medium ISP#2

ISP #2 (BD 277010)..... 38 g
DI Water..... 1000 mL

Autoclave at 121°C for 15 minutes.

Note: Broth Media must be made from scratch

To make media from scratch use formulation below

Malt Extract 10.0 g
Yeast Extract..... 4.0 g
Dextrose..... 4.0 g
*Agar..... 20.0 g
DI Water..... 1000 mL

*Omit agar for broth formulation

Adjust for pH 7.3 +/- 0.2

Autoclave at 121°C for 15 minutes.