

**ATCC medium: 2003 Czapek yeast extract agar**

K<sub>2</sub>HPO<sub>4</sub> .....1.0 g  
Czapek Concentrate (see below)....10.0 ml  
Yeast extract.....5.0 g  
Sucrose.....30.0 g  
Agar.....15.0 g  
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.

*Czapek Concentrate:*

NaNO<sub>3</sub> .....30.0 g  
KCl.....5.0 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....5.0 g  
FeSO<sub>4</sub> . 7H<sub>2</sub>O .....0.1 g  
ZnSO<sub>4</sub> . 7H<sub>2</sub>O .....0.1 g  
CuSO<sub>4</sub> . 5H<sub>2</sub>O .....0.05 g  
Distilled water.....1.0 L