

ATCC medium: 2031 LB + 1.5 mg/ml trimethoprim

LB Medium (ATCC Medium 1065, see below) plus 1.5 mg/ml trimethoprim

Calculate the required amount of trimethoprim reagent, dissolve in solvent, and filter-sterilize.

Add to sterile basal medium.

ATCC Medium 1065:

Tryptone (BD 211705).....	10.0 g
Yeast Extract (BD 212750).....	5.0 g
NaCl.....	10.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.