

ATCC medium: 2075 Medium "A" for Alkaliphiles

Glucose.....	10.0 g
Peptone, Bacto (BD 211677)...	5.0 g
Yeast extract.....	5.0 g
KH ₂ PO ₄	1.0 g
MgSO ₄ . 7H ₂ O	0.2 g
NaCl.....	40.0 g
Na ₂ CO ₃	10.0 g
Agar.....	20.0 g
Distilled water.....	1.0 L

Adjust medium, if necessary, for final pH 9.6. Autoclave at 121C for 15 minutes.