

ATCC medium: 2102 Marine medium with 0.5% alginate

Note: Use either Marine Agar 2216 (BD 212185) or Marine Broth 2216 (BD 279110). Add Sodium alginate.

Marine Broth (BD 279110).....37.4 g
Sodium alginate.....5.0 g
Distilled water.....1.0 L

Add sodium alginate very, very slowly to dissolve. Then boil, dispense, and autoclave at 121C for 15 minutes.

Marine Agar (BD 212185).....55.1 g
Sodium alginate.....5.0 g
Distilled water.....1.0 L

Add sodium alginate very, very slowly to dissolve. Then boil, dispense, and autoclave at 121C for 15 minutes.