

**ATCC medium: 2162 Horikoshi-I medium**

Glucose.....10.0 g  
Polypeptone.....5.0 g  
Yeast extract.....5.0 g  
K<sub>2</sub>HPO<sub>4</sub> .....1.0 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....0.2 g  
Agar.....15.0 g  
Distilled water.....900.0 ml

Autoclave at 121C for 15 minutes.

After autoclaving, aseptically add 100.0 ml of sterile 10% Na<sub>2</sub>CO<sub>3</sub> to the medium.

Adjust for final pH of 10.0.