

ATCC medium: 2168 *Halomonas* medium with elevated NaCl

NaCl.....	250.0 g
Casamino Acids (BD 223050).....	7.5 g
Proteose Peptone No. 3 (BD 211693).....	5.0 g
Yeast extract.....	1.0 g
Sodium citrate.....	3.0 g
MgSO ₄ . 7H ₂ O	20.0 g
K ₂ HPO ₄	0.5 g
Fe(NH ₄) ₂ (SO ₄) ₂ . 6H ₂ O	0.05 g
Distilled water.....	1.0 L

Adjust pH to 7.0 with KOH.
Autoclave at 121C for 15 minutes.